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# CareMat® with replaceable cable

The pressure-sensitive bed mat  
for dementia patients and persons at risk of falling.

## Translation of the original instructions

### 1 Safety instructions

These instructions must be completely read through prior to the commissioning and use of the CareMat®. They should be kept safely for future use. Warranty claims will be invalid in the event of any damage to the product, to goods or persons that is due to failure to observe these instructions.

#### Instructions for nursing staff

- The CareMat® operates purely as a support for the nursing staff, but is not a substitute for them.
- The functional correctness of the CareMat® is to be checked each time that it is used.

#### Notes for the installer


- The electrical installation work and commissioning of the product may only be carried out by a qualified technician
- The product must not be connected to mains voltage (see "Different Types of Electrical Connections").

### 2 Product description

#### 2.2 Areas of application and function

The CareMat® was primarily designed for use in old people's and care homes, hospitals and psychiatric wards as well as for sheltered residential accommodation. The necessary care and supervision of patients with dementia and of those who are at risk of falling represents a significant increase in care expenditure for care practitioners. The CareMat® is the ideal support: Risks of falling and the dangers which arise from wandering or loss of orientation can be significantly reduced.


The CareMat®, as a pressure-sensitive sensor, is installed directly in front of the care-patient or sick person's bed and is connected to the internal call system. As soon as the person leaves the bed and steps onto the CareMat® this switches on and sends out an alarm signal to the call system via the RCL07.


CareMat A01  
1100 x 700mm

- ① Replaceable cable for connecting to the call system
- ② Open cable end
- ③ Housing with terminals or universal circuit-board
- ④ Chamfered, wheelchair-accessible edges
- ⑤ Non-slip surface with stud texture design
- ⑥ Carrying handle


### 3 Day to day operation

Position the CareMat® centrally in front of the bed so that avoiding the mat is impossible. Route the cable under the bed in order to avoid the risk of tripping over it. Insert the plug into the correct plug - socket device.

If the CareMat® should not be needed, it can be pushed beneath the bed without the cable connection to the call system having to be disconnected.

Remove, CareMat® before moving the bed around. By doing this, damaging the mat or the cable will be avoided.

### 4 Cleaning and disinfecting

Disconnect the connection to the call system prior to any cleaning.

- CareMat® can be cleaned using commercially available disinfecting and cleaning agents with max. 70 vol% alcohol. No corrosive agents should be used.
- CareMat® is neither autoclavable nor can it be disinfected using steam.

### 5 Transport and storage

- During transport any significant bending of the CareMat® should be avoided.
- Store the CareMat® in a dry, clean and flat location.
- Do not place any heavy loads on the CareMat® and do not stack more than 20 packed mats on top of each other.

## 6 Different electrical connections - with terminals circuit-board

The CareMat A01C (C=cable) includes a terminal circuit-board in the housing.

The correct plug for connecting to the alarm system and its terminal assignment are to be obtained from the relative system manufacturer.

## 7 Different electrical connections - with universal circuit-board



The CareMat A01P (P= circuit-board) includes a universal circuit-board in the housing.

The universal circuit-board can be adapted to your needs by altering the switch positions. The wiring on the universal circuit-board must not be removed.

The connections are protected against wrong polarity of input voltage.

### Variable resistance R

The universal circuit-board has two connector sockets for commercially-available wired resistors. A call resistor that is specifically for the call system for line monitoring can be installed in this. The connecting wiring should have a diameter of between 0.38 and 0.56 mm. The socket spacing should be 14 mm.

⚠ After the electrical installation, the installer must perform a functional test.

## 8 Maintenance and fault finding

The CareMat® is maintenance-free. If the cable should however be pulled out, it should then simply be replaced: open housing, replace and connect up cable. The cable ties around the connecting cable must be stripped using a cable tie tool to avoid any slipping. Close the cover and lightly tighten the screws.

If despite this the function is not being provided, the following checks should be made:

- Is the call system operating and is the CareMat® correctly connected?
- Is the plug, the plug connection or the cable damaged?

If the points mentioned above have been checked out, but the functions of the CareMat® are, despite this, still not being provided, you should consult a qualified expert.

## 9 Disposal

The CareMat® can be disposed of in the household refuse.

## 10 Technical specifications

<b>Article</b>	CareMat A, K / B	<b>Cable</b>	Length 3 m, Ø 5.1 +/- 2 mm, cable cross-section ≤ 0.14 mm <sup>2</sup>
<b>Dimensions (LxWxH)</b>	1100x700x9 mm / 700x400x9 mm, height with housing 15 mm	<b>Type of cable</b>	4-Pol
<b>Weight [kg]</b>	8.2, 6.9 /3.2	<b>Signal output</b>	Make contact or resistor connection
<b>Material</b>	Polyurethane (PU)	<b>Using the circuit-board</b>	0 mA at 10 mA
<b>Colour</b>	Mat: care yellow, Cable: white	<b>Electrical load capacity</b>	Contact voltage: max. 25 mA at 24 V
<b>Surface</b>	Stud texture design	<b>Type of protection</b>	IP 54 without cable end / connecting plug
<b>Weight of persons</b>	At least 10 kg	<b>Ambient temperature</b>	0 up to +55° C
<b>Chemical resistance</b>	Water (room temperature), body fluids, disinfectant, max. 70 vol% alcohol	<b>Conformity to standards</b>	The CareMat® with terminal board is just subject to the General Product Safety Directive 2001/95/EC and does not therefore require CE labelling. The CareMat® with circuit board comes under the EMC Directive 2004/108/EC and therefore has CE conformity labelling.

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