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REF 2154 - 2156

NORMATIVE

This product is a Class I medical device and complies with Regulation (EU) 2017/745 of the European Parliament and of the Council of April 5, 2017 on Medical Devices.

This product should not be disposed of in the waste container. There are specific collection points where you can deposit them. For the disposal of the container and the product, strictly comply with the legal regulations of your community.

INTENDED USE

- Support product indicated to be used by people who have functional deficits in one or both lower extremities caused by pathologies such as arthritis and osteoarthritis, neurological diseases or some type of trauma or surgery.
- Indicated in patients with: broken bones, postoperative, osteoarthritis, plaster or the need to provide greater balance, both statically and dynamically.

CORRECT USE

- It must be held in the hand opposite the affected leg, either due to weakness or surgery. Two crutches can be used for greater stability.
- The elbow should be slightly flexed when resting on the crutches.
- The handle should be at the level of the wrist and the greater trochanter (hip bone).
- Never hang bags or other objects on the crutches.
- Do not make any modifications (mechanical alterations) to the crutches so as not to compromise its safety of use.
- Crutches should always be used for their intended purpose. They must not be used for any other purpose. For example, they should not be used as a lever to lift an object, such as push-up bars, etc.
- The user must not exceed the maximum weight capacity of the product (indicated on the product label)

PRECAUTIONS

- Before each use, you must make sure that the product is in perfect conditions.
- Before each use, make sure the tip is in perfect condition. If it shows wear, contact the healthcare facility/store where you purchased the product immediately for replacement or repair.
- Do not use the product in case of visible defect. In the event of missing, wrong fitting or damaged parts; seek advice from the healthcare facility/store where you purchased the product.
- Increase caution when using on carpets, wet surfaces, slippery surfaces (leaves, gravel, sand or snow).
- Do not use on frozen surfaces.
- Keep the product away from heat sources as it could cause deterioration or flammability. Do not expose the product to direct heat sources such as stoves, heaters, radiators, etc.

REF	Product reference		Manufacturing Date
MD	Medical Device	LOT	Lot number
CE	CE Certificate		Instructions for use
	Manufacturer		Precaution
	Distributor		Importer

If instructions, warnings, and notes in the user manual are not followed, it may cause physical injury to you or others, or damage to the product.

CONTRAINDICATIONS

If the product is not adjusted to the correct height, it can create bad posture vices that over time can cause injuries and muscle contractures in the back, arms and hands.

ASSEMBLY AND INSTALLATION OF ADJUSTABLE CRUTCHES

To adjust the height: As this cane has a pin, first open the anti-kickback locking tab. Then remove the locking clip located on the side of the tube by pressing the outer edges of the clip with your thumbs. Choose the desired height by extending or shrinking the tube. When you reach the desired height, place the pin in the selected hole. After adjusting the crutch, do not forget to close the anti-eject lock tab again. Remember to set both English poles at the same height.

CONSERVATION AND WASHING RECOMMENDATIONS

Clean and dry periodically. Remove dirt with a soft, damp cloth and neutral soap. Dry later. Do not use aggressive or abrasive chemical cleaners; they can damage the material of the product. Also, do not use greasy products that can degrade materials or make it slippery.

When not in use, store the product in a dry place at room temperature.

SHELF LIFE

The useful life is influenced by its conditions of use and storage and by the regularity of its maintenance. By imperatively respecting the data provided in this instruction manual, it is possible to keep the product in perfect condition for 24 months.

WARRANTY

SUCESORES DE MANUEL GARCIA S.L. guarantees all its products, provided that they have not been manipulated or altered in their original configuration. It does not guarantee those products whose characteristics are altered due to misuse, deficiencies or breakage of any kind. The product is guaranteed for 24 months from the date of purchase for manufacturing defects. Excluded from the guarantee are damages caused by improper use, accident or natural wear of parts such as: ferrules and other plastic components, footrests, upholstery, brakes, wheelchair tires, and similar parts, depending on the type of medical device purchased.

WARNING

- If you observe any deficiency, anomaly or failure, immediately notify the vending establishment to proceed with its change or repair.
- Any serious incident related to the product must be reported to Sucesores de Manuel García, S.L. and to the competent authority of the Member State in which the user and/or the patient are established.

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