**DE****Schlafbrille – Blau-Blocker – UV530**

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2024.01.31

Read this package leaflet carefully before you use the Sleep Glasses.

- Keep this leaflet. You may need to read it again.
- Consult your retailer, MelaMedic or your doctor if you want to know more.
- Speak to your doctor if you experience any side effects.

Content of the leaflet

1. Effect and application
2. What you need to know before you use Sleep Glasses
3. How to use Sleep Glasses
4. Side effects
5. Storage and cleaning
6. Safety
7. Verification
8. Explanation of Symbols
9. Additional Information

1. EFFECT AND USE

The sleep glasses are intended for the prevention and relief of sleep disorders. They are also used to correct disrupted sleep patterns, for example, in relation to night shifts and jet lag.

The use of sleep glasses stimulates the secretion of the melatonin, a neurotransmitter from the pineal gland in the brain, which promotes sleep, sleep quality, and other physiological processes.

2. WHAT YOU NEED TO KNOW BEFORE YOU USE SLEEP GLASSES

GB

Activities

You can easily watch TV, read, work at your computer and participate in most of your usual activities in the evening, whilst you are wearing your Sleep Glasses.

Adaptation time and regularity

It can take up to a week to get used to your new sleep glasses and achieve the desired effect. This is due to the fact that you are in the process of resetting your body's inner clock. In some cases, the adaptation time can be up to 30 days and then it is a matter of being patient and following the instructions carefully.

It is important to put the Sleep Glasses on at the same time every evening in order to maintain the same setting of your internal clock.

For example, if you wear the Sleep Glasses regularly every evening from 21, then the melatonin cycle will typically stop at 07 the next morning. If you wake up too early and feel refreshed and rested, then you can try to wear the Sleep Glasses from 22.00 instead.

If you follow the instructions, you will – like the vast majority – already experience a positive effect on your sleeping pattern within the space of three – four days.

There are a few people whose bodies are only able to produce little or no melatonin. For these people the effect of using Sleep Glasses can be limited.

3. HOW TO USE SLEEP GLASSES

Sleep disorders

Sleep disorders can be caused by imbalances in the natural circadian rhythm, as the natural production of melatonin is primarily inhibited by blue light in the light spectrum and partially by green light in the light spectrum (UV-550). Blue light is often referred to as light pollution (artificial lighting, TV and screens).

Sleep Glasses (UV-530) block all the blue light, and by wearing them 1-3 hours before bedtime, the melatonin cycle can be advanced so that the secretion of melatonin has decreased or stopped before waking up the next morning. It is easier to wake up and feel refreshed when the melatonin level is low.

If you wear the Sleep Glasses as directed 1-3 hours before bedtime and stay in the dark whilst you are sleeping, then you can achieve a cumulative melatonin production of the normal 8-10 hours.

For example, if you wear the Sleep Glasses regularly every evening from 21, then your melatonin cycle will typically end the next morning at 07.

When you retire to your darkened bedroom to sleep, you must switch off the light before you take the Sleep Glasses off. If you have to get out of bed during the night e.g. to go to the toilet, then you must either do this in the dark or else put on the Sleep Glasses whilst doing so.

If you expose your eyes to normal white light (which contains blue light) even briefly, then you risk resetting the production of melatonin and you will lose the whole effect that you have

built up with the Sleep Glasses over the course of the evening.

GB

Night shift

A combination of Sleep Glasses and light therapy can help night workers to regulate their circadian rhythms and sleep patterns so that they manage to sleep well during the day and can feel refreshed and energetic for their night work.

If you work from e.g. from 23 to 07, it is recommended to use light therapy before starting work. It will be easier to feel fresh and alert as the light from light therapy lamps emits extra blue light, which inhibits melatonin production. It is also beneficial to add light therapy while working.

As soon as you leave your night shift, you must put on the Sleep Glasses, which block all the blue light, which will lead to an increased production of melatonin, stimulating sleep. If, e.g. you sleep from 09 to 17, then you should put on the Sleep Glasses as soon as you leave your night shift. The Sleep Glasses get the body to slow down and get it ready for sleep.

When you wake up at 17, you must immediately perform light therapy for about 30 minutes. And then start all over again.

Jet lag

Most people that travel over long distances experience and are plagued by jet lag.

The symptoms of jet lag could be fatigue, indigestion, insomnia, sleep disorders, concentration difficulties, bad moods, low energy and general malaise.

Jet lag is caused by a disruption of the body's inner clock, which is a little collection of

braincells that control the circadian rhythm and a series of biological functions. This includes mealtimes and sleep. The body's inner clock is constructed according to a regular rhythm between daylight and darkness. The inner clock goes out of sync when daylight and darkness are at completely different times in the new time zone to those the body is used to. That is why you get jet lag.

The more time zones you pass through, especially when flying eastwards, the worse the experience of jet lag is certain to be. Jet lag happens when you travel too far and too fast. It is due to the body's inner clock being adjusted to your usual time zone.

By using the Sleep Glasses on the day of departure and possibly the following day, according to a simple time schedule, you can reduce the symptoms of jet lag significantly.

Based on a few specific details, you can get a free Jet Lag Guide on www.melamedic.com, with instructions on how you should respectively seek or avoid light. In this way, you will be able to adjust to a new time zone within 1-2 days instead of the normal 5-6 days.

4. SIDE EFFECTS

There are no reported side effects to using Sleep Glasses when they are used according to the instructions for use. When you change your circadian rhythm with the aid of the Sleep Glasses then it can be expected that some people may experience dizziness in the first days.

If you experience any side effects, please inform MelaMedic immediately.

5. STORAGE AND CLEANING

GB

The glasses must be stored in the accompanying spectacle case. This applies to both storage and transport. Avoid pressure on the case and climatic heat when transporting.

The glasses are best cleaned in warm water to remove grease, dirt, sand, cosmetics, etc. Polish the glasses afterwards with a soft, clean cloth made out of cotton or microfibre.

6. SAFETY

⚠ Do not use the Sleep Glasses when driving a vehicle or operating machinery as they may pose a security risk due to incorrect colours from lamps and signs in addition to increased fatigue.

⚠ It is crucial for the user to follow the instructions to achieve the desired effect and avoid unintended incidents.

⚠ In severe cases of illness, the Sleep Glasses should only be used in consultation with your doctor.

⚠ Sleep Glasses with any visible damage to the frames should not be used due to the risk of injuries to the eyes and face.

⚠ The spectacle lenses are made out of polycarbonate and the eyeglass frame is made out of nylon. This means that the glasses should not be used close to open flame or hot objects, because they are flammable and can melt. Avoid climatic heat.

⚠ The Sleep Glasses should last a long time. It depends on how they are used and stored. The glasses should be changed if the lenses are scratched or dull due to improper storage, use or cleaning.

7. VERIFICATION

The Taiwan Spectacles Industry Association continuously monitors the glasses using light spectrum analysis. This ensures that the glasses meet the requirement for filtering out light up to 530 nanometers.

This product has been tested and manufactured in accordance with the requirements of Medical Device Regulations 2017/745/EU for Class 1 medical devices.

8. EXPLANATION OF SYMBOLS

This mark indicates that the product is manufactured by the manufacturer specified by the symbol.

This mark indicates that the product complies with all EU regulations and directives.

This symbol indicates the importance of reading the package leaflet before using the medical device.

This mark provides cautionary information for the user of the medical device.

9. ADDITIONAL INFORMATION

GB

The product identification is located on the inside of the frame:

melamedic

www.melamedic.com TYPE: UV530

Medical Device Class I - 2017/745/EU
For therapeutical use only
MM-XXX / YYMMDD

Manufacturer and holder of the marketing authorization for this product:

MelaMedic ApS
Trindsovej 6
8000 Aarhus
Denmark
Reg. No.: DK40641947

The Declaration of Conformity is available through:

www.melamedic.com

This package leaflet was last revised on:
2024.01.31

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2024.01.31

All rights reserved

All of the content in this leaflet including product images, trademarks, illustrations, design, icons, photographs, artwork and images are copyrights, trademarks and other intellectual property owned and controlled by MelaMedic.

References**Bipolar Disorder, 2016 May;18(3):221-32**

Blue-Blocking glasses as additive treatment for mania: a randomized placebo-controlled trial.

J Adolesc Health, 2015 Jan;56(1):113-9

Blue Blocker glasses as a countermeasure for alerting effects of evening light-emitting diode screen exposure in male teenagers.

Chrono Physiology and Therapy, 2013 January 4

Treatment of attention deficit hyperactivity disorder insomnia with blue wavelength light-blocking glasses.

Med Hypotheses, 2008;70(2):224-9

Dark therapy for bipolar disorder using amber lenses for blue light blockade.

Journal of Pineal Research, 2006 Aug;41(1):73-8

Blue Blocker glasses impede the capacity of bright light to suppress melatonin production.

Chronobiology International, 2006;23(1-2):475-84

Thermoregulatory effects of melatonin in relation to sleepiness.

Journal of Clinical Endocrinology & Metabolism, 2005 May;90(5):2755-61

Blocking low-wavelength light prevents nocturnal melatonin suppression with no adverse effect on performance during simulated shift work.

Journal of Neuroscience, 2001 Aug 15;21(16):6405-12

Action spectrum for melatonin regulation in humans: evidence for a novel circadian photoreceptor.

Journal of Physiology, 2001 Aug 15;535(Pt 1):261-7

An action spectrum for melatonin suppression: evidence for a novel non-rod, non-cone photoreceptor system in humans.

Hormone Research, 1998;49(3-4):118-24

Effect of seasonal changes in daylength on human neuroendocrine function.

melamedic®
C I R C A D I A N E Y E W E A R