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Run programmed sequence
Run a custom sequence to measure your improvement

Hand/eye coordination

Measure your reaction time and count hits within a given time

Program Sequences

Create and modify sequences to challenge yourself and your friends

Split time test

Set the FitLights on a line and measure the time between your passes

Export training data

Export results from your training sessions and analyse them in Fitlight or 3rd part tools

PROGRAM SEQUENCE

1-8

Created: 05.28.2014 - 02:04:22
Modified: 05.28.2014 - 03:37:46
Sequence steps: 8
Lights required: 8

Fitlight_sequence_04_08_2014_063404

Created: 04.08.2014 - 06:34:12
Modified: 06.21.2014 - 12:52:12
Sequence steps: 3
Lights required: 2

Fitlight_sequence_07_02_2014_080610

Created: 07.02.2014 - 08:06:12
Modified: 07.02.2014 - 08:06:57
Sequence steps: 5
Lights required: 4

attack

Created: 07.02.2014 - 05:00:18
Modified: 07.02.2014 - 05:25:20
Sequence steps: 29
Lights required: 5


3. Vælg en eksisterende sekvens at redigere, eller tryk på "New" for at lave en ny sekvens

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
PERFORMANCE & TRAINING

- Random training**
Adjust fitLight settings and let the application generate randomized sequences
- Run programmed sequence**
Run a custom sequence to measure performance
- Hand/eye coordination**
Measure your reaction time and count hits within a given time
- Export training data**
Export results from your training sessions and analyse them in Fitlight or 3rd part tools

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PROGRAMMED SEQUENCE SETUP

Running Sequence:
Name: attack
Created: 07.02.2014 - 05:00:18
Modified: 07.02.2014 - 05:25:20
Sequence steps: 29
Lamps required: 5
Maximum execution time: ~46s

+ REPEATS **+** GLOBAL DELAY
- 2 **-** 0.00 s

Start by light Save data
 Singular Random start

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