
Robin™

Robin Mover™

Invacare® Robin™

Robin™ , Robin Mover™

en	Ceiling Hoist Service Manual	3
da	Løftemotor Servicemanual	33
fr	Lève-personne sur rail Manuel de maintenance	65
it	Sollevatore a soffitto Manuale per la manutenzione	99
nl	Plafondlift Servicehandleiding	133
no	Takmontert løfter Servicehåndbok	165
sv	Taklyft Servicemanual	197

DEALER: Keep this manual.
The procedures in this manual MUST be performed by a qualified technician.

Yes, you can.®

© 2015 Invacare® Corporation

All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

Hilti is a registered trademark of Hilti Corporation.

fischer is a registered trademark of fischer fixings UK Ltd.

Loctite is a registered trademark of Henkel Corporation.

SPAX is a registered trademark of SPAX International GmbH & Co.

Contents

1 General	4
1.1 About this manual	4
1.2 Safety information	4
2 Components	5
2.1 Main parts of the ceiling hoist	5
2.2 Accessories	5
2.3 Spare parts	6
3 Service	7
3.1 Installing the ceiling hoist to the rail system	7
3.1.1 Installing Robin	7
3.1.2 Installing Robin Mover	8
3.2 Dismounting the ceiling hoist	9
3.2.1 Dismounting Robin™	9
3.2.2 Dismounting Robin Mover™	10
3.3 Replacing the hand control	12
3.4 Removing the hoist cover	12
3.5 Replacing Strap and Hook Assemblies	12
3.6 Replacing the battery	13
3.7 Replacing the PCB (Printed Circuit Board)	14
3.8 Replacing the fuse	15
3.9 Adjusting the solenoid until the hoist traverses	15
3.10 Adjusting Loose Strap Switch	16
3.11 Replacing the Twin Speed Switch	18
3.12 Replace Contact Disc for Robin Mover™	18
4 Maintenance	20
4.1 Charging the Battery	20
4.2 Battery Maintenance	20
4.2.1 Storage Characteristics	21
4.2.2 Storage Temperature	21

4.2.3 Storage Time	21
4.2.4 Storage Humidity	21
4.2.5 Capacity Recovery	22
4.3 LOLER Statement	22
4.4 Safety Inspection Checklist	23
5 Troubleshooting	24
5.1 Troubleshooting Table	24
5.2 Indicator light and audio signal	26
5.3 Checking Loose Strap Function	26
5.4 Checking Twin Speed Functionality	26
5.5 Taking Voltage Readings	27
6 Technical data	29
6.1 The rail system	29
6.2 Dimensions and weight	29
6.3 Electrical System	30
6.4 Environmental conditions	30
6.5 Electromagnetic compliance (EMC) information	31
6.5.1 Guidance and manufacturer's declaration	31

I General

I.1 About this manual

This service manual contains information about assembly, adjustment and advanced maintenance of the product. In order to ensure safety when handling the product, read the manual carefully and follow the safety instructions.

For pre-sale and user information, please see the user manual. Find the user manual on Invacare's website or contact your Invacare representative (see addresses at the end of this manual).

Symbols in this manual

In this manual, hazards are indicated by symbols. The symbols are accompanied by a signal word that indicates the severity of the risk.

WARNING

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

CAUTION

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.

IMPORTANT

Indicates a hazardous situation that could result in damage to property if it is not avoided.

Tips and Recommendations

Gives useful tips, recommendations and information for efficient, trouble-free use.

Gives information about the required tools.

This product complies with Directive 93/42/EEC concerning medical devices. The launch date of this product is stated in the CE declaration of conformity.

I.2 Safety information

WARNING!

The procedures in this service manual, must be performed by a specialized dealer or qualified service technician.


- Do not handle this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manuals, installation manuals or instruction sheets supplied with this product or optional equipment.

The information contained in this document is subject to change without notice.

2 Components

2.1 Main parts of the ceiling hoist

(A)	Connection for hand control
(B)	Control button UP
(C)	Directional arrows on Robin Mover™
(D)	Control button DOWN
(E)	Connection for external back-up battery
(F)	Indicator light
(G)	Lifting straps
(H)	Red emergency strap
(I)	Strap hooks

(J)	Retainer for hand control (accessory)
(K)	Hand control

2.2 Accessories

CAUTION!

Compatibility of slings and hanger bars / strap hooks

Invacare® uses a "Loop and Coat Hanger Bar / Strap Hook System" as do many other manufacturers. Therefore other suitable patient transfer systems (slings), manufactured by other companies, can be used on the Invacare patient lift range as well.

However we do recommend:

- A risk assessment is always to be carried out by a professional prior to issuing lifting equipment. It is important that the Task, Individual, Load, Environment and Equipment are considered in the risk assessment.
- Always choose the sling design and size according to the patient's weight, size and physical ability whilst considering the type of transfer to be carried out.
- Do only use slings that are suitable for a "Loop and Coat Hanger Bar / Strap Hook System".
- Do not use slings for "Keyhole Hanger Bar" or for "Tilting Frame Hanger Bar" designs.

Available accessories

- 2- and 4-point hanger bar including flexiscope version
- Stretcher for horizontal lift
- Scale including hanger bar
- Holder for hand control
- 24 V charger

- Charger set (24 V charger and holder for hand control)
- Extension straps and extension cables for hand control
- External back-up battery
- Charger for external back-up battery
- Sling models for “Loop and Coat Hanger Bar System”:
 - Full body support slings – without head support
 - Full body support slings – with head support
 - Slings for dress/toileting – with or without head support
 - Slings for amputee
 - Slings for walking training

2.3 Spare parts

An electronic spare parts catalogue is available on Invacare’s website. If you don’t have access, contact your Invacare representative (see addresses at the end of this manual) for a printed spare parts catalogue.

3 Service

3.1 Installing the ceiling hoist to the rail system

Robin™ requires a trolley, guiding the ceiling hoist in the rail system. Robin Mover™ has the trolley built-in and can be installed directly to the rail system.

3.1.1 Installing Robin

Tools: 4 mm Allen key; 17 mm socket wrench

1. To install the trolley, remove the end stop that is attached to the rail system.

In some countries/ local areas, there is a demand for double securing of the end stops. If an extra securing screw is installed, remove this screw before.


2. Insert the trolley into the rail (no matter which way, as the trolley is symmetric).
3. Re-insert the end stop.


In some countries/ local areas, there is a demand for double securing of the end stops. If an extra securing screw is required, install this screw. See also End Stop.

4. Locate the hoist trolley at the Ø 32 mm round recess in the rail.


5. Lift the hoist and push the suspension pin through the trolley opening.

6. Turn the hoist by 90 degrees to fix it in the trolley.
7. Connect the hand control.
8. Adjust the length of the emergency stop/lowering cord to within the user's reach.

If there is excessive resistance, or if the trolley unintentionally rolls by itself due to imbalance in the rail suspensions, adjust the friction brake using a 4 mm Allen key.


3.1.2 Installing Robin Mover

The trolley is built-in on the hoist, and can therefore not be pre installed in the rail.

Tools: 4 mm Allen key; 17 mm socket wrench

1. Remove the end stop.

2. Insert the wheels of the hoist into the rail.



 The direction the hoist faces does not matter during installation.

3. Install the end stop.

 Tighten the bolt to 30 Nm.

4. If too much resistance is experienced when the hoist is pulled, or the hoist moves too freely and unintentionally, adjust the friction brake.



5. Connect the hand control.
6. Adjust the length of the emergency stop/lowering cord to within the user's reach.

3.2 Dismounting the ceiling hoist

3.2.1 Dismounting Robin™

Remove the ceiling hoist from the trolley

1. Lift the hoist slightly and push the suspension pin out of its holder in the trolley.
2. Turn the hoist by 90 degrees and pull it out of the trolley.

Remove the trolley from the rail

1. Remove the end stop.
2. Slide the trolley out of the rail.

3.2.2 Dismounting Robin Mover™

1. Pull the emergency strap to prevent accidental activation of the ceiling hoist after it has been dismantled from the rail.
2. Remove the end stop.

3. Slide the wheels of the hoist out of the rail.

4. Install the end stop.


 Tighten the bolt to 30 Nm.


5. If removing the hoist from the rail motor, perform the following steps:

 When replacing the battery or performing service, it is also necessary to remove the hoist from the rail motor.

- a. Remove the screws securing the plastic shells from the rail motor.

- b. Pull the plastic shells away from each other to remove them from the rail motor.

- c. Remove the M6 screw securing the rail motor to the hoist.

- d. Separate the rail motor from the hoist.

CAUTION!


Risk of damage

It is possible to remove the hoist from the rail motor without removing the rail motor from the rail system.

To avoid damage to the hoist:

- Always support the hoist when removing the M6 screw.

3.3 Replacing the hand control


1. Slide and hold locking lever **A** away from the aperture **C**.
2. Pull the plug **B** out of the aperture.
3. Insert plug of the new hand control into aperture. Make sure, the plug is in the correct orientation with the slot upwards (see image).
4. Release lever **A**, so that the gate closes and hold the plug in the aperture.

3.4 Removing the hoist cover

Before removing the hoist cover, pull the red emergency strap to turn off the hoist and remove the hoist from the track system. Refer to:

- 3.2 Dismounting the ceiling hoist, page 9

 Tools: Screw driver, slotted

1. Remove protection caps **A**. Use fingers or, if required, use a screw driver carefully.
2. Untighten screws **B**, using a screw driver.
3. Remove the cover **C** from the hoist.


To close the hoist cover, tighten the screws and attach the protection caps again.

3.5 Replacing Strap and Hook Assemblies

 Tools: 10 mm socket key; 5 mm hex/Allen key

1. With the ceiling hoist in the track, lower the straps to their lowest point using the hand control.
2. Pull the red emergency strap to perform an emergency stop. Do not perform emergency lowering.
3. Remove hoist from track. Refer to 3.2 Dismounting the ceiling hoist, page 9.
4. Remove the cover. Refer to 3.4 Removing the hoist cover, page 12.

5. Reset the emergency stop by pressing the reset button.
6. Check that the strap **A** is completely unwound as shown in the image. If not, wind the straps out using the UP/DOWN button on the hand control.
7. Align the strap retaining bolts **B** with the aperture in the chassis using the UP/DOWN button on the hand control.

8. Pull the red emergency strap again.
9. Remove strap bolts, using a 10 mm socket key and a 5 mm hex/Allen key.
10. Remove the strap through guide cylinder.
11. Insert a new strap through guide cylinder and assemble using the strap retaining bolts (Torque: 4 Nm). Make sure that the new strap is fed over the slack strap switch wire.
12. Re-mount the cover, re-assemble the hoist to the track and reset the emergency stop.

3.6 Replacing the battery


CAUTION!

Static electricity can damage the PCBs inside the hoist.

Before working in the inside of the hoist:

- Pull the red emergency strap to turn off the hoist.
- Remove the hoist from the track system. (Refer to 3.2 Dismounting the ceiling hoist, page 9)
- Remove the cover from the hoist. (Refer to 3.4 Removing the hoist cover, page 12)
- An anti-static mat must be used when working inside the ceiling hoist.

1. Unplug the wires **A**.
2. Remove the battery pack **B**.
3. Insert a new battery pack.

Use only Invacare recommended batteries.

4. Plug in the wires (A).
5. Re-assemble the hoist and reset the emergency stop.


3.7 Replacing the PCB (Printed Circuit Board)


CAUTION!

Static electricity can damage the PCBs inside the hoist. Before working in the inside of the hoist:

- Pull the red emergency strap to turn off the hoist.
- Remove the hoist from the track system. (Refer to 3.2 Dismounting the ceiling hoist, page 9)
- Remove the cover from the hoist. (Refer to 3.4 Removing the hoist cover, page 12)
- An anti-static mat must be used when working inside the ceiling hoist.

1. Remove the battery. Refer to steps 1 and 2 in chapter 3.6 Replacing the battery, page 13.
2. Remove all plugs (B), (C), (D), (E) and (F) from the PCB (A).
3. Release the PCB (A) from the retaining holders by gently pressing down the locking clips and levering the PCB slightly away from the chassis plate.

4. Lift the PCB up to remove.
5. Place the new PCB into the chassis and press it slightly against the chassis plate until the retaining holders snap into place.
6. Connect the plugs (B), (C), (D), (E) and (F) to the new PCB.
7. Re-insert battery, re-assemble the hoist and reset the emergency stop.

3.8 Replacing the fuse


CAUTION!

Static electricity can damage the PCBs inside the hoist.

Before working in the inside of the hoist:

- Pull the red emergency strap to turn off the hoist.
- Remove the hoist from the track system. (Refer to 3.2 Dismounting the ceiling hoist, page 9)
- Remove the cover from the hoist. (Refer to 3.4 Removing the hoist cover, page 12)
- An anti-static mat must be used when working inside the ceiling hoist.

1. Remove the hand control. Refer to chapter 3.3 Replacing the hand control, page 12
2. Remove the battery. Refer to steps 1 to 3 in chapter 3.6 Replacing the battery, page 13.
3. Gently pull the fuse **A** straight out away from the holder **B**.

4. Insert new fuse.

Make sure only a fuse with correct rating 10A 250V slow blow is used.

5. Re-insert battery, re-assemble the hoist, insert the hand control and reset the emergency stop.

1586747-A



3.9 Adjusting the solenoid until the hoist traverses

Robin Mover™ only.


4 mm Allen key; 7 mm spanner

1. Pull the red emergency strap.
2. With the ceiling hoist remaining on the track, remove both plastic shells:

- a. Remove the screws securing the plastic shells from the rail motor.
- b. Pull the plastic shells away from each other to remove them from the rail motor.

- Using a cloth, protect the solenoid plunger whilst clamping with grips at the same time.
Using a 7 mm spanner raise the nut (A) by turning anti-clockwise 2 complete turns.

- For testing, restrain the hoist with one hand whilst pressing the Left or Right button on the hand control.
There should be good contact between the white drive wheel and the underside of the track. If not, repeat step 3.
- Reassemble the plastic shells. Ensure the correct orientation, with clearance for the plate above the solenoid to move up and down freely. See images in step 2.

3.10 Adjusting Loose Strap Switch

CAUTION!

Static electricity can damage the PCBs inside the hoist.


Before working in the inside of the hoist:

- Pull the red emergency strap to turn off the hoist.
- Remove the hoist from the track system. (Refer to 3.2 Dismounting the ceiling hoist, page 9)
- Remove the cover from the hoist. (Refer to 3.4 Removing the hoist cover, page 12)
- An anti-static mat must be used when working inside the ceiling hoist.


2 mm Allen key

1. Place the hoist on a surface (such as the corner of a table) so that the hoist is supported yet the straps and hand control are free hanging. Ensure the hoist is safe and cannot slip/fall of the table.
2. Locate the loose strap switch **A** and visually check if there is clearance between collet **C** and chassis.

- a. If not: whilst retaining the geometry of the slack strap wire **B** and collet arm, loosen the collet grub screw **D** with a 2 mm Allen key.
 - b. Create clearance between collet and chassis.
 - c. Retighten the collet grub screw.
3. Test the loose strap function. Refer to chapter 5.3 Checking Loose Strap Function, page 26. This test can be carried out with the hoist remaining on the table. If the fault is rectified, re-assemble the hoist. If not, continue with step 4.

4. Adjust tension on the slack strap wire:

- a. Loosen the collet grub screw **D** with a 2 mm Allen key. Note the spring **F** will now pull the collet arm **C** against the switch so it is open circuit.
 - b. Lift the strap **E** off the slack strap wire **B** and lift up the slack strap wire.
 - c. Retighten the collet grub screw **D**.
5. Test the loose strap function. Refer to chapter 5.3 Checking Loose Strap Function, page 26. This test can be carried out with the hoist remaining on the table.
 6. If the fault is rectified, re-assemble the hoist and reset the emergency stop.

3.11 Replacing the Twin Speed Switch

CAUTION!

Static electricity can damage the PCBs inside the hoist.


Before working in the inside of the hoist:

- Pull the red emergency strap to turn off the hoist.
- Remove the hoist from the track system. (Refer to 3.2 Dismounting the ceiling hoist, page 9)
- Remove the cover from the hoist. (Refer to 3.4 Removing the hoist cover, page 12)
- An anti-static mat must be used when working inside the ceiling hoist.

T10 Torx screw driver

1. Place the hoist on a surface (such as the corner of a table) so that the hoist is supported yet the straps and hand control are free hanging. Ensure the hoist is safe and cannot slip/fall off the table.
2. Remove the two yellow cables (A) from the twin speed switch (B).

3. Remove the two screws (C).
4. Attach the new twin speed switch with screws (C). Connect the yellow cables (A). Ensure, they are connected to the bottom and top terminals.
5. Re-assemble the hoist and reset the emergency stop.

3.12 Replace Contact Disc for Robin Mover™

CAUTION!

Static electricity can damage the PCBs inside the hoist.


Before working in the inside of the hoist:

- Pull the red emergency strap to turn off the hoist.
- Remove the hoist from the track system. (Refer to 3.2 Dismounting the ceiling hoist, page 9)
- Remove the cover from the hoist. (Refer to 3.4 Removing the hoist cover, page 12)
- An anti-static mat must be used when working inside the ceiling hoist.


Screwdriver

1. Remove the cable tie **(A)**.
2. Remove the contact disc cables **(B)** from the Mover PCB, using a terminal screwdriver. Note brown cable top and blue cable bottom.
Feed the cables through the chassis plate.

3. Carefully separate and remove the complete contact disc **(C)** from the suspension base **(D)** by releasing the self-adhesive tape.
4. Remove, clean off the remains of the self adhesive tape.

5. Fit the new self adhesive tape **(E)** and the new contact disc **(C)**.

6. Feed the cables **(B)** through the chassis plate and connect them to the Mover PCB.
7. Re-assemble the hoist and reset the emergency stop.

4 Maintenance

4.1 Charging the Battery

If the charge status becomes low, an audible indicator sounds and the indicator lamp on the bottom of the ceiling hoist flashes (see also 5.2 Indicator light and audio signal, page 26). When this happens, there is usually sufficient power to lower the patient.

It is not possible to use the hoist while the battery is charging.

1. Finish the current lift in progress
If the current lift in process discharges the battery completely, an external backup battery may be used to temporarily power the lift. See Using the external back-up battery.
2. Move the ceiling hoist to the charging station.
3. Place the hand control into the charging station.
A yellow indicator lamp on the charger indicates the hand control is connected. The indicator lamp illuminates red while charging.
4. When fully charged, the indicator lamp on the charger turns green.
As long as the ceiling hoist is not in use, leave the hand control in the charging station.
5. To use the ceiling hoist, remove the hand control from the charger.

Keeping the batteries charged guarantees hoist functionality and maintains the battery to ensure a long lifetime.

If the audible indicator beeps and the blue indicator on the hoist flashes after the hoist has been fully charged and only been in use for a short period of time, the battery may need exercise or replacement. See 4.2 Battery Maintenance, page 20.

If the ceiling hoist will not be used for more than 4 weeks, Invacare recommends disconnecting the charger and pulling the emergency stop.

The charger is connected to the mains with a mains plug. Make sure that the mains plug is accessible and can be unplugged if required.

4.2 Battery Maintenance

If the audible indicator beeps and the blue indicator on the hoist flashes after the hoist has been fully charged and only been in use for a short period of time, the battery may need to be exercised.

To enable the battery to charge fully, the battery should be exercised from time to time.

Following a long period of storage, it could also be necessary to exercise the battery.

It is also suggested that the battery is exercised in situations where the hoist is used infrequently, and/or where the hand control is placed in the charger for long periods of time. In this situation, there will always be capacity enough for a few lifts, but we recommend to exercise the battery as soon as possible afterwards.

Exercise of the battery can be done by charging and discharging the hoist completely several times.

1. Charge the battery until the LED turns green.
2. Remove the hand control from the charger.
3. Raise and lower the straps with a load for one minute.
4. Wait ten minutes.
5. Repeat Steps 3–4 until the battery is completely discharged.


Continue this process after the audible alert sounds and the blue light flashes, indicating the battery is low.

6. Repeat Steps 1–5 two to three times.

4.2.1 Storage Characteristics

Essentially all rechargeable battery cells gradually discharge over time whether they are used or not. This capacity loss is typically due to slow parasitic reactions occurring within the cell. As such, the loss rate (self-discharge rate) is a function of the cell chemistry and the temperature environment experienced by the cell. Due to the temperature sensitivity of the self-discharge reaction, relatively small differences in storage temperature may result in large discharges which are difficult or impossible to reverse. Cell and battery storage issues of concern to most application designers relate either to the speed with which the cells lose their capacity after being charged or the ability of the cells to charge and discharge “normally” after storage for some period of time. In both situations, general guidelines developed for nickel-cadmium cells work acceptably for nickel-metal hydride cells.

4.2.2 Storage Temperature

As already mentioned, the self-discharge reaction rate increases with higher temperatures. Prolonged storage of the battery material deteriorating faster; leakage performance will also deteriorate, resulting in a reduced battery lifetime. It is recommended that, for long storage, batteries should be kept at room temperature or below (0-30°C or 32–86°F).

4.2.3 Storage Time

As the battery loses energy during storage, the voltage also drops. In general, the battery capacity loss due to self-discharge during storage can be recovered by recharging. If the battery is stored for over six months, it is advisable to cycle the battery several times to resume the battery capacity. Use good inventory practices (first in, first out) to reduce time cells spend in storage.

4.2.4 Storage Humidity

Leakage and rusting of metal parts are accelerated in high humidity environments, especially those with correspondingly high

temperatures. The recommended humidity level for battery storage is a maximum of 60% RH.

4.2.5 Capacity Recovery

In normal practice, stored cells will provide full capacity on the first discharge after removal from storage and charging with standard methods. Cells stored for an extended period or at elevated temperatures may require more than one cycle to attain pre-storage capacities. Consultation with the manufacturer is recommended if prolonged storage and rapid restoration of capacity is requested.

4.3 LOLER Statement

The UK Health and Safety Executive's Lifting Operations and Lifting Equipment Regulations 1998, require any equipment that is used in the workplace to lift a load be subject to safety inspection on a six monthly basis. Please refer to the HSE web site for guidance www.hse.gov.uk.

The person responsible for the equipment must ensure adherence to LOLER regulations.

4.4 Safety Inspection Checklist

IMPORTANT!

– Perform this inspection every 12 months.

Rail system

- Check that all brackets are tightened correctly to the ceiling/wall. Retighten, if necessary.
- Check reinforcement. If not possible, do a load test with 1.5 x rated load for 20 minutes.
- Check and retighten end stops.
- Systems with track connections: Check if the connection has become misaligned during use. Re-align if needed.
- Check trolley and traverse trolley.
- Check trolley wheels for smooth movement and abnormal noise.
- Check friction break adjustment on trolleys. Adjust if necessary.
- Do a load test with rated load through the entire system.
- Place a sticker with a new date for the next inspection.

Rail system with transit coupling connections:

- Check the lock fixtures:
 - secure to tracking
 - each claw moves independently
 - correctly aligned with tracking slots
 - no obstructions to prevent correct function
 - safety plate is fitted.
- Check traverse trolley mounting fixtures.
- Check tracking alignment. Re-align if necessary.
- Function test transit coupling / lock fixture (10 times - disable left / right claw x 5).
- Check that gap between track and transit coupling is max. 3 mm.
- Function test under load (5 cycles with Safe Working Load).
- Check inlay plate is fitted (L and M rail).

Ceiling hoist

- Check suspension and cross pin.
- Check steel and plastic gearwheels for wear and damage. If gearwheels are dry, apply gearwheel grease.
- Check loose strap sensor and adjust if necessary.
- Check straps for wear and damage.
- Check guide cylinder and end stops.
- Check the lifting hooks, screws and pins for wear and damage.
- Check that a sticker with maximum load is placed on the hoist.
- Do a full lifting cycle with rated load — must be effortless and without jarring sound.
- Place a sticker with a new date for next inspection.

Electrical connections

- Check that all cables are intact.
- Check for rupture on cables and on isolation.
- Check micro switches for clamping and mechanical function.
- Check the charger for damage / function.
- Check hand control for damage and function. Make sure that the hand control is connected correctly.
- Check the emergency stop.
- Check the emergency lowering.

Serial Number:	
Inspection Date:	
Inspected by:	
Comments:	

5 Troubleshooting

5.1 Troubleshooting Table

SYMPTOMS	PROBLEM	SOLUTION
Ceiling hoist does not respond to the hand control's buttons and up/down buttons on the ceiling hoist are not responding.	Emergency stop activated	Reset the emergency stop. (Refer to user manual)
	System does not have power or battery is discharged	Charge the battery or use the emergency battery. Refer to 4.1 Charging the Battery, page 20.
	Fuse is blown	Replace the fuse. Refer to 3.8 Replacing the fuse, page 15.
Ceiling hoist does not respond to the hand control's buttons, but up/down buttons on the ceiling hoist are responding.	Hand control not connected properly	Check that the hand control cable is properly plugged into the hoist. Refer to 3.3 Replacing the hand control, page 12.
	Hand control defective	Replace the hand control. Refer to 3.3 Replacing the hand control, page 12.
Lifting straps will only move down and not up.	Battery nearly discharged	Charge the battery or use the emergency battery. Refer to 4.1 Charging the Battery, page 20.
	Load is too high	Reduce the load. The hoist is equipped with an over-current safety system.
Lifting straps will only move up and not down.	Hand control defective	Replace the hand control. Refer to 3.3 Replacing the hand control, page 12.
	The loose strap function activated	Check the loose strap function. Refer to 5.3 Checking Loose Strap Function, page 26.
Ceiling hoist speed is the same with and without load.	Non-load quick movement function activated permanently	Check the quick movement switch. Refer to 5.4 Checking Twin Speed Functionality, page 26.
Lifting straps are twisted and do not straighten out.	Lifting hooks do not turn freely enough	Clean the hooks or replace the straps.
	Straps are frayed	Replace the straps. Refer to 3.5 Replacing Strap and Hook Assemblies, page 12.

SYMPTOMS	PROBLEM	SOLUTION
Ceiling hoist does not traverse the rails.	Rails require service or cleaning	Inspect the rail system
	Driving wheels worn or damaged (Robin Mover™ only)	Adjust driving wheels.
	Hand control is defective	Replace the hand control.
Ceiling hoist emits excessive noise when activated.	Bearings, gearwheels or motor defective	Replace the hoist.
Ceiling hoist will not charge	System has no power	Place the hand control in the charging station and recharge the ceiling hoist
		Check that the power supply is switched on and connected to the charging station.
No indicator light on hoist	System does not have power or battery is discharged.	Check that the ceiling hoist has power supply and that the battery is recharged.
		Check that the power supply is switched on and connected to the charging station.
		Place the hand control in the charging station and recharge the ceiling hoist.
The emergency strap does not turn hoist off	Hoist requires service.	Contact an Invacare representative for repair.
The emergency strap does not lower the hoist.	Hoist requires service.	Contact an Invacare representative for repair.

If problems are not remedied by the suggested means, contact an Invacare representative.

5.2 Indicator light and audio signal

Hoist

Status	Indicator light	Audio signal
Off / Stand by Hoist traversing (Robin Mover™ only)	Off	No
Hoist in use (Up/Down)	Blue	No
Low battery	Blue Flashing	Yes

The blue indicator light on the underside of the hoist turns off automatically after 30 seconds.

Wall Mounted Charger

Status	Indicator light	Duration
Battery not connected	Yellow	n/a
Initializing and analysis of battery	Yellow	12–14 seconds
Fast charge	Red	5 min — 3 hours
Charged (trickle charge)	Green	n/a
Fault condition (battery defect)	Alternating red and green	n/a

5.3 Checking Loose Strap Function

A correctly adjusted loose strap switch ensures that the straps stop lowering when they touch an obstruction.



1. With the ceiling hoist remaining on the track, check there are no obstructions under either strap hook.
2. Place an obstruction under one of the strap hooks, then press and hold the down button on the hand control, Expected result: When the strap hook touches the obstruction, the strap stops lowering.
3. If the strap does NOT stop lowering, adjust the loose strap switch. Refer to chapter 3.10 Adjusting Loose Strap Switch, page 16.
4. Repeat steps 2 and 3 for the second strap.

5.4 Checking Twin Speed Functionality

1. With the hoist on the track and both straps free from obstructions, press and hold the up button with no load applied. Note speed (40 mm/s) and sound of motor.
2. Whilst continuing to hold the up button, apply a load of at least 10 kg (e.g. by pulling down on both straps). Note the speed reduction (30 mm/s) and motor sound change.
3. If there is no speed reduction, replace the twin speed switch. Refer to chapter 3.11 Replacing the Twin Speed Switch, page 18.

5.5 Taking Voltage Readings

<p>Read Feeding Voltage</p>	
<p>Take battery voltage feeding PCB. Expected result: 23 to 29 V</p>	
<p>Take voltage feeding main switch (two middle pins) Expected result: 23 to 29 V</p>	
<p>Take voltage feeding the lift motor whilst pressing the UP or DOWN button on the hand control. Expected result: 23 to 29 V</p>	

<p>Read Feeding Voltage</p>	
<p>Take voltage feeding the hoist mover PCB. (Robin Mover™ only) Expected result: 23 to 29 V</p>	
<p>Take voltage feeding the copper rings whilst pressing the LEFT or RIGHT button on the hand control. (Robin Mover™ only) Expected result: 23 to 29 V</p>	
<p>Take voltage feeding the trolley mover PCB whilst pressing the LEFT or RIGHT button on the hand control. (Robin Mover™ only) Expected result: 23 to 29 V</p>	

Read Feeding Voltage	
<p>Take voltage feeding the trolley mover motor whilst pressing the LEFT or RIGHT button on the hand control. (Robin Mover™ only)</p> <p>Expected result: 23 to 29 V</p>	
<p>Take voltage feeding the solenoid whilst pressing the LEFT or RIGHT button on the hand control. (Robin Mover™ only)</p> <p>Expected result: 23 to 29 V</p>	

6 Technical data

6.1 The rail system

All rail profiles are manufactured in extruded aluminium.

The rails are available with white powder coating.

The rails are available in three different heights; S= 68 mm, M=114 mm and L=183 mm, providing a free span of 2 m, 4 m and 8 m, respectively.

The profile width of all three profiles is 58 mm. The top of the profiles are shaped to permit the coupling of the rails to the installation brackets.

On delivery, the rail for the hoist will feature a milled groove on its bottom. This groove is used for hoist installation and is called keyhole.

Possible installations

Single rail can be installed on the ceiling or onto the wall and floor. The rail position in the room can be either parallel to a wall or angular in relation to a wall.

The traverse system can be installed on the ceiling or onto the wall and floor. The free span between the two parallel rails can be up to 8000 mm.


Curves are available in the S-profile and can be installed on the ceiling. They are available in angles of 30°, 45°, 60° and 90°.

Special solutions

Installation with pendling brackets is used in the following cases: When the height of the ceiling varies in the lengthwise direction of the rail, or when the rail is to be countersunk.

Angular installation is used when the rail has to be installed in other angles than 90°.

6.2 Dimensions and weight

Dimensions

	Robin™	Robin Mover™
Hoist length (A)	470 mm	470 mm
Hoist width (B)	250 mm	250 mm
Hoist height (C)	203 mm	266 mm
Maximum lifting range (D)	2.5 m	2.5 m

Weights

	Robin™	Robin Mover™
Maximum lifting capacity (Safe Working Load)	200 kg	200 kg
Total weight without sling	7.3 kg	10 kg

6.3 Electrical System

	Robin™	Robin Mover™
Maximum current input	0.9 A	
Protection class ¹	Entire device: IPx4 (ceiling hoist: IP24; hand control: IP44)	
Insulation class	Class II equipment	
	Type B applied part Applied part complying with the specified requirements for protection against electrical shock according to IEC60601-1.	
Sound level	50 – 55 dB (A)	
Number of lifts per charge	65 lifts of 0.5 m with 80 kg	
Intermittent (periodic motor operation)	10%, max. 2 min. / 18 min.	
Battery	NiMH Battery 24V / 2.0 Ah	
Charging time	Up to 3 hours at 20°C	
Lifting speed with 200 kg load	2.1 cm/s	2.1 cm/s
Traverse speed in the tracks	n/a	14 cm/s

Charger

Voltage output charger	24 V DC, max. 250 VA
Voltage supply charger	230 V AC, 50 Hz
Protection class ¹ charger	IP41 or IP67 (depending on version)

¹ See product label and label on each electric device for correct protection class. The lowest IP classification decides the overall classification of the device.

- IPx4: Protected against water splashed from any direction.
- IP24: Protected against objects larger than 12.5 mm and protected against water splashed from any direction.
- IP44: Protected against objects larger than 1 mm and protected against water splashed from any direction.
- IP41: Protected against objects larger than 1 mm and protected against vertically dripping water.
- IP67: Dust-tight and protected against temporary immersion into water.

6.4 Environmental conditions

	Storage and transportation	Operation
Temperature	-30 °C to +50 °C	+5 °C to +40 °C
Relative humidity	10 % to 70 %	20 % to 90 % at 30 °C, not condensing
Atmospheric pressure	700 hPa to 1060 hPa	

6.5 Electromagnetic compliance (EMC) information

Medical Electrical Equipment needs to be installed and used according to the EMC information in this manual.

This equipment has been tested and found to comply with EMC limits specified by IEC/EN 60601-1-2 for Class B equipment.

Portable and mobile RF communications equipment can affect the operation of this equipment.

Other devices may experience interference from even the low levels of electromagnetic emissions permitted by the above standard. To determine if the emission from the lift is causing the interference, run and stop running the lift. If the interference with the other device operation stops, then the lift is causing the interference. In such rare cases, interference may be reduced or corrected by the following:

- Reposition, relocate, or increase the separation between the devices.

6.5.1 Guidance and manufacturer's declaration

Test	Basic Standard	Compliance Level / Frequency Range
Mains Port Conducted Emissions	EN 55011:2007	Group I Class B
Radiated Emissions	EN 55011:2007	Group I Class B
Mains Harmonics	EN 61000-3-2:2006 + A1:2009 + A2:2009	Class A

Test	Basic Standard	Compliance Level / Frequency Range
Voltage Fluctuations	EN 61000-3-3:2008	dc%/dmax/d(t)/Pst
Electrostatic Discharge	EN 61000-4-2:1995	±2/4/6kV Contact, ±2/4/8kV Air
Radiated Field Immunity	EN 61000-4-3:2006	3V/m, 80% 1kHz AM, 80MHz - 2.5 GHz
EFT/Bursts	EN 61000-4-4:2004	±2kV L/N/E/LNE
Surge Immunity	EN 61000-4-5:2006	±0.5/1kV L-N, ±0.5/1/2kV L-E/N-E
Conducted RF Immunity	EN 61000-4-6:2007	3V, 80% 1kHz AM, 0.15MHz - 80MHz
Magnetic Field Immunity	EN 61000-4-8:1993	3A/m, 50Hz & 60Hz
Dips & Interruptions	EN 61000-4-11:2004	-95%/-60%/-30%

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Invacare Ireland Ltd,
Invacare AB
Invacare UK Operations Limited