# Lisa

|   |     |
|---|-----|
| <b>DE</b> Gebrauchsanweisung .....        | 3   |
| <b>EN</b> Instructions for Use.....       | 29  |
| <b>FR</b> Instructions d'utilisation..... | 55  |
| <b>IT</b> Istruzioni per l'uso .....      | 81  |
| <b>ES</b> Instrucciones de uso.....       | 107 |
| <b>PT</b> Manual de utilização .....      | 133 |
| <b>NL</b> Gebruiksaanwijzing .....        | 159 |
| <b>DA</b> Brugsanvisning .....            | 185 |
| <b>PL</b> Instrukcja użytkowania .....    | 211 |
| <b>HU</b> Használati utasítás .....       | 237 |
| <b>CS</b> Návod k použití.....            | 263 |
| <b>HR</b> Upute za uporabu .....          | 287 |
| <b>BS</b> Uputstvo za upotrebu .....      | 311 |
| <b>RU</b> Руководство по применению ..... | 335 |

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## Instructions for Use for LISA Rehab Folding Buggy/ Stroller (models from 2005, HR3213.../HR3214...)

| Contents   | Page      |
|--|-----------|
| <b>1 General Information</b> .....                                   | <b>30</b> |
| <b>1.1 Preface</b> .....   | <b>30</b> |
| <b>1.2 Intended Use</b> .....  | <b>30</b> |
| <b>1.3 Liability</b> .....   | <b>30</b> |
| <b>1.4 CE conformity</b> .....                                       | <b>31</b> |
| <b>1.5 Service and Repairs</b> .....                                 | <b>31</b> |
| <b>2 Safety Instructions</b> .....                                   | <b>31</b> |
| <b>2.1 Warning Symbols</b> .....                                     | <b>31</b> |
| <b>2.2 General Safety Instructions</b> .....                         | <b>32</b> |
| <b>2.3 Warning Symbols and Type Plates</b> .....                     | <b>34</b> |
| <b>3 Delivery and Preparing the LISA for Use</b> .....               | <b>34</b> |
| <b>4 Adjustment and Adaptation Possibilities</b> .....               | <b>37</b> |
| <b>5 Accessories</b> .....   | <b>41</b> |
| <b>5.1 Abductor</b> .....  | <b>41</b> |
| <b>5.2 Storage Bag</b> .....   | <b>42</b> |
| <b>5.3 Armrests</b> .....  | <b>42</b> |
| <b>5.4 Lap belt</b> .....  | <b>42</b> |
| <b>5.5 Buggyboard</b> .....  | <b>43</b> |
| <b>5.6 Canopy</b> .....  | <b>43</b> |
| <b>5.7 Five-point belt</b> .....                                     | <b>44</b> |
| <b>5.8 Grabrail and grabrail padding</b> .....                       | <b>45</b> |
| <b>5.9 Hook and loop foot straps with instep pads</b> .....          | <b>46</b> |
| <b>5.10 Headrest</b> .....   | <b>46</b> |
| <b>5.11 Swivel lock</b> .....  | <b>46</b> |
| <b>5.12 Padding</b> .....  | <b>47</b> |
| <b>5.13 Rain cape</b> .....  | <b>48</b> |
| <b>5.14 Slipsacks</b> .....  | <b>48</b> |
| <b>5.15 Lateral guides with variable seat width reduction</b> .....  | <b>49</b> |
| <b>5.16 Groin strap</b> .....  | <b>49</b> |
| <b>6 Using your LISA for transportation in a motor vehicle</b> ..... | <b>51</b> |
| <b>7 Maintenance and Cleaning</b> .....                              | <b>51</b> |
| <b>8 Technical Data</b> .....  | <b>52</b> |
| <b>9 Re-use / Disposal</b> .....                                     | <b>53</b> |
| <b>9.1 Information on Re-use</b> .....                               | <b>53</b> |
| <b>9.2 Information on Disposal</b> .....                             | <b>53</b> |

# 1 General Information

## **INFORMATION**

Date of the last update: 2012-08-15

- Please read this document carefully.
- Follow the safety instructions.

### **1.1 Preface**

With the LISA rehab folding buggy/ stroller you have decided for an innovative and functional state-of-the-art aid.

This manual shows how you can quickly and easily make use of all functions. The instructions on necessary maintenance and cleaning shall help you enjoy this product for a long time to come.

Should you have further questions, please contact your authorized dealer.

The design, as described in these Instructions for Use, is subject to technical alterations without notice.

### **1.2 Intended Use**

The prescription of rehab buggies/ strollers is generally indicated for the severely disabled, especially for children and young people who are unable to walk, need to be transported and are unable to operate a wheelchair themselves.

Inability to walk or very serious walking impediment due to:

- Paralysis
- Loss of limbs
- Defective and/or deformed limbs
- Joint contractures
- Joint defects
- Other diseases

### **1.3 Liability**

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according to the instructions for use.

The manufacturer is not responsible for damages caused by components and spare parts not approved by the manufacturer. Repairs must be carried out exclusively by authorized dealers or by the manufacturer.

## 1.4 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

## 1.5 Service and Repairs

Service and repairs on the LISA may only be carried out by Ottobock approved dealers. Should any problems arise, please contact the supplier of your LISA.

Any necessary repairs will be made exclusively with authentic Ottobock spare parts there.

### Note!

For repair and service, the following tools are required:

Allen wrench, sizes: 3 - 6 mm

Wrench, sizes: 10 mm, 13 mm, 19 mm

Your authorized Ottobock dealer:

## 2 Safety Instructions

### 2.1 Warning Symbols

#### Danger!

Warning messages regarding possible risks of accident or injury.

#### Attention!

Warning messages regarding possible technical damages.

#### Note!

Information regarding operation of the device.

#### Note!

Information for service staff.

#### Attention!

Please read the Instructions for Use first!

## 2.2 General Safety Instructions

**Danger!**

We would like to point out that using the rehab folding buggy/ stroller beyond the normal use, may be dangerous. The product is unsuitable for jogging, running, skating etc. Swivelling front wheels may begin to wobble at higher speeds and lead to an abrupt stop, causing the LISA to tip over. Therefore, please only push the buggy/ stroller at walking speed. It is in no case permissible to let go of the push handles while pushing the buggy/ stroller or to thrust it away from you.

**Danger!**

Use caution near flammables or fire, especially when holding lit cigarettes. The seat and backrest padding might catch fire

**Danger!**

Static stability is at 15° inclination. Attaching heavy bags or other weight to the push handles can adversely affect stability.

**Danger!**

Always strap the child in the rehab folding buggy/stroller.

**Danger!**

The product may be used only on solid level ground.

**Danger!**

Before leaving the LISA and before getting into and out of it, always lock the brakes.

**Danger!**

Stairs should be ascended/descended only with assistance from attendants. If ramps or elevators are available, they should be used. If such devices are missing two helpers must carry the wheelchair over the obstacle.

**Danger!**

Pay particular attention when on slopes and inclines to prevent

- the child from falling out of the LISA
- the LISA from tipping over
- the LISA from rolling away

**Danger!**

Park the buggy only on level surfaces. If you have to park the LISA on an incline, make sure that you bring the backrest to an upright position. There is a risk that the LISA can tip over backwards when in a reclined position.

**Danger!**

When your child reaches for objects in front, to the side or behind the LISA, make sure that it does not lean out too far, since the shift in the center of gravity might cause the LISA to tilt or tip over.

**Danger!**

Please keep packaging material away from children. Plastic packaging presents the danger of suffocation.

**Danger!**

After you have finished all adjustments, snugly re-tighten the screws that have been loosened.

**Danger!**

Never leave your child unattended in the LISA, even if the belts are in use and the brakes are engaged.

**Danger!****Using your Ottobock Product as a Seat for Transportation in Motor Vehicles for Disabled Passengers**

We recommend that, wherever and whenever possible, users of a wheelchair, mobility base with seating shell or buggy transfer to the seats installed in the motor vehicle and use the corresponding vehicle restraint systems, because this is the only way to ensure optimal protection of the passengers in case of an accident.

Your Ottobock product Lisa is permissible for transporting passengers in motor vehicles when using Ottobock safety components as well as appropriate restraint systems. For more information, please refer to our instructions for use manual "**Using your wheelchair/ mobility base with seating shell or buggy/stroller as a seat for transportation in motor vehicles**", order no. 646D158.

**Attention!**

The correct function of the rehab folding buggy/ stroller must be checked before each use.

**Attention!**

Use your rehab buggy/stroller properly. For instance, do not drive into obstacles (including steps, curbs) without braking.

**Attention!**

The footplate must not be used as a stepstool to get in or out of the buggy/stroller.

**Attention!**

Tire pressures have a big effect on the ease of handling. Correctly inflated tires considerably improve the maneuverability of your LISA.

The air pressure should be at least 2 bar (200 kPa).

**Attention!**

Observe the maximum load (body weight including any additional load) of the LISA.

Model 3213 = 50 kg (110 lbs) / model 3214 = 60 kg (132 lbs)

**Attention!**

For clearing obstacles such as steps and curbs, tilt the LISA onto the rear wheels (to ascend, pull the Lisa backwards; to descend, slowly lower the Lisa forward).

**Attention!**

Lift the rehab folding buggy/stroller only by components which are firmly attached (by welding or screw connection).

-Front frame tube above the front wheels

-Push handles

**Note!**

The rehab folding buggy/ stroller has been approved for the transport of only one child.

**Note!**

All accessories and add-on components will reduce the remaining maximum load capacity

**Attention!**

Please read the Instructions for Use first!

Make sure you are familiar with the product and how to operate it before the first use. You are responsible for the safety of your child. You may compromise the safety of your child, if you do not observe these instructions. However, it is not possible to avoid any unforeseen circumstances and situations.

Common sense, caution and circumspection are factors that cannot be provided by the product but are expected from the person using the buggy/ stroller. The user of the buggy/ stroller and its accessories should understand all instructions. The operation of the LISA should be explained to any person who uses it and its accessories.

If these Instructions for Use are unclear, need to be explained or if you have further questions, please contact an Ottobock authorized dealer.

### 2.3 Warning Symbols and Type Plates

| Label / Nameplate | Explanation  |
|-------------------|--|
|                   | <ul style="list-style-type: none"> <li>A Type designation</li> <li>B Manufacturer's article number</li> <li>C Maximum load capacity (see Section „Technical Data“)</li> <li>D Manufacturer / address / country of manufacture</li> <li>E Serial number</li> <li>F Manufacturing date</li> <li>G European Article Number / International Article Number</li> <li>H Read the Instructions for Use prior to using the product.</li> <li>I CE marking – product safety in conformity with the EU Directives</li> </ul> |
|                   | <p>Anchoring point/securing point for securing the rehab device in motor vehicles for the transportation of disabled persons.</p>  |

### 3 Delivery and Preparing the LISA for Use

**Danger! – Risk of injury!**

Do not hold on to moveable parts other than the ones described, especially the folding scissor construction.

**Attention!**

Damage caused by objects lying upon the appliance. If appliances are transported in the folded state, the weight of objects lying upon the appliance can cause deformations and result in problems when unfolding the appliance. Never put heavy objects on a folded rehabilitation device.

In the original package you will find the following components:

- LISA in folded position with mounted footrest holders
- Footrest for LISA (not mounted)
- Instructions for use and tools required
- Accessories as ordered

First, carefully remove the packaging and place the buggy/ stroller in front of you (Fig. 1).

**To unfold the LISA please proceed as follows:**

- First disengage the transport lock (Fig. 2).
- Hold on to a push handle to lift the LISA (Fig. 3).
- Place one foot behind a rear wheel as illustrated (Fig. 3).

- Tip the LISA slightly onto a rear wheel (Fig. 3).
- Push the folding lever (ball grip) down in semi-circular motion and unfold the LISA (Fig. 4).
- Make sure the slides rest properly on the stops

- Press down the folding lever at the backrest (ball grip) until it reaches the stop (Fig. 5).

**When assembling the footrest for the first time, follow the instructions under section 4 "Adjusting the Footrest".**

**In everyday use, proceed as follows:**

- Fold down the footrest and slide the footrest receiver slightly upwards. Then turn the footrest bar outwards to inwards (Fig. 6).

- Join the footplate and the remaining footrest bar together (Fig. 7).
- Put the footplate's notch onto the holder on the footrest bar (Fig. 8).

- Engage the lock with a quarter turn of the locking knob so that the bolt locks into the footrest (Fig. 9).

**Danger!**

After you have finished all adjustments, snugly re-tighten the screws that have been loosened.

**Attention!**

Make sure the footplate has securely locked by giving it a few short tugs.

**To fold the LISA, you proceed as follows:**

- Remove the accessories from the LISA.
- Unlock the footrest in the reverse order as described and fold it away to the side.
- Release the brake.
- Bring the backrest to upright position and pull the lock open.
- Proceed in reverse order as described under "unfolding"
- Finally, secure with the transport lock.

## 4 Adjustment and Adaptation Possibilities

### Footrest

Upon delivery, the footrest and the two footrest bars are not mounted.

To carry out the assembly, proceed as follows:

- Connect the free footrest bar to the footplate (Fig. 8+9).
- For assembly, turn the footrest receivers at the front frame to the front and insert the footrest/ footrest bars into the footrest receivers from below (Fig. 10).
- Before tightening the screws, make sure that the footrest bars are at the same height and the footrest receivers are locked to prevent distortion.

### Setting the Lower Leg Length


After loosening the clamping screws on the left and right side of the footrest receiver, the footrest can be adjusted in height (Fig. 10).

Starting from a lower leg length of 34 cm (13"), the footrest receiver at the anterior frame must be set deeper.

Loosen and remove the holding screw (Fig. 11) and slide the footrest receiver to the second bore hole. Retighten the screws at the footrest receiver. Repeat these modification steps at the opposite side of the front frame.

### Setting the Footrest deeper

The depth of the footrest can be adjusted to two positions (Fig. 12).

### Setting the Seat Depth

First loosen and remove the screws (Fig. 13, pos. A) on both sides of the seat bottom. Afterwards, carry out assembly at the desired seat depth (Fig. 14).

#### **Danger!**

After you have finished all adjustments, snugly re-tighten the screws that have been loosened.

### Application of the Brake

- To apply the brake, press the brake pedal (Fig. 15) down with your foot.
- To release the brake (Fig. 16), pull the brake pedal upwards with your foot.

### Backrest Inclination

- First lock the brake.
- Secure the backrest with one hand (Fig. 17).
- Pull the release cord a few centimeters (Fig. 17) to release the (Fig. 18) angle adjustment of the backrest.

- Adjust the desired backrest angle while keeping the release cord pulled (Fig. 19 + 20). Let go of the release cord to let the bars on the left and right side lock in place (Fig. 18).

### **Cover of the Folding Scissor Construction**

The folding scissor construction of the backrest is covered with a pocket-type cover (Fig. 21). Only put small and soft objects into the pocket, because at maximum back inclination these parts may become very close to the frame tubes.

### **Tip-Assist**

The tip-assist may be screwed onto the left side or the right side of the posterior frame. To tilt the rehab folding buggy/ stroller, place one foot on the tip-assist (Fig. 22) and pull the push handles in your direction.

## 5 Accessories

### Note!

All accessories not installed by the manufacturer must be installed by trained technicians.

The following notes on installation are listed here for technician use.

Before installation of many accessories, the leather strips must be cut open. To do this first mark the leather. Next pre-punch the start and end of the cut with a hole punch and finally make a cut that connects two holes.

Straps on accessories are usually extra long to accommodate most positions.

Excess strap ends can be cut off, and the cut edge melted together with a flame (cigarette lighter) in order to prevent fraying.

### Note!

The accessories which always have to be removed before folding the rehab buggy/ stroller are marked with a dot (⊗).

### 5.1 Abductor

The abductor is screwed to the folding scissor construction under the seat (Fig. 23).

To ease getting into and out of the LISA, fold the abductor forward (Fig. 24).

To unlock, fold the abductor in direction of the seat bottom, press it down and keep it pressed down while folding the abductor forward.

### Attention!

The storage bag can be loaded up to a maximum of 5 kg (11 lbs).  
The storage bag must be emptied before folding the LISA.

### 5.2 Storage Bag (Fig. 25)

The storage bag is attached to the frame under the seat and folds together with the rehab buggy/stroller.

### 5.3 Armrests (Fig. 26)

Screw the armrests onto the left and right side of the front frame. They can be mounted at two heights. The armrests must be flipped up before folding the buggy/stroller.

### 5.4 Lap belt

The standard equipment of the LISA includes a simple universal lap belt (Fig. 27), which can be cut to length, if required.

Alternatively, lap belts are offered as an accessory in different sizes (Fig. 28).

The lap belt runs over the front of the backrest upholstery and through the two eyelets on the lower edge of the backrest. The lap belt is opened by the quick-release buckle. Depending on the model, the length of the lap belt is adjusted at the quick-release buckle, or the two-piece lap belt is adjusted at the hook and loop closure. The strap ends must overlap by at least 15 cm (6").

### 5.5 (⊗) *Buggyboard* (Fig. 29)

The buggyboard is intended for transporting a second child behind the folding buggy/stroller. It is provided with two spring-mounted rollers. The buggyboard is attached to the back frame. It can be flipped up or detached using the quick closures. Please observe the Instructions for Use/Assembly Instructions enclosed with the product (Accessory not offered in USA).

#### **Danger!**

The buggyboard can be loaded with up to a maximum of 20 kg and is unsuitable for scooter riding.

### 5.6 (⊗) *Canopy*

The canopy (Fig. 30) is clipped to the upper ends of the push handles. Pull the canopy over the backrest and close the hook and loop closures around the sides of the push handles.



If it rains, pull the enclosed transparent raincover (Fig. 31) over the canopy and down over the footrest assembly. The push handles remain uncovered.

When using the raincover, only an upright backrest position is possible.

## 5.7 (⊗) **Five-point belt**

The five-point belt can be used with or without padding. It consists of the shoulder belts, the lap belt and the centre belt. All of the belts meet in the closure buckle. The buckle can only be closed if the belt tabs are inserted in the correct order (Fig. 32).

When using the five-point belt without padding (Fig. 33a), the shoulder straps are attached with one strap that is run around the backrest (Fig. 33b).

When using the five-point belt together with the padding (Fig. 34a), proceed as follows:

- 1) At the upper edge of the back rest, open the hook-and-loop closure on the padding (Fig. 34b, Pos. 1). Push the padding ahead slightly.
- 2) Open the zippers on the top and bottom of the padding (Fig. 34b, Pos. 2) and pass the shoulder belts through the gap between the zippers from the front.
- 3) Set the desired shoulder height by adjusting the zipper position.
- 4) Pass the five-point belt behind the padding and around the back rest; fasten the closure buckle (Fig. 33a).
- 5) At the top edge of the back rest, fasten the hook-and-loop closure on the padding (Fig. 34b, Pos. 1).

The **lap belt** runs over the front of the backrest upholstery and through the two eyelets (Fig. 35) at the lower end of the backrest. The length of the lap belt is adjusted via the belt buckle.

The **center strap** prevents the lap belt from slipping upwards. The center strap is attached to the corresponding loop at the seat bottom (Fig. 36) and adjusted in length via the buckle.

**Note!**

The center strap and lap belt must not be removed before folding the rehab folding buggy/stroller.

### 5.8 (⊗) **Grabrail and grabrail padding**

To mount the grabrail, the receivers must be screwed to the seat tubes (Fig. 37).

Insert both ends of the grabrail into the receivers until the stop. The spring buttons must be pressed down.

The arch of the grabrail must point towards the seat.

To remove the grabrail (Fig. 38), press down the spring buttons and keep them pressed.

Afterwards, pull the grabrail out of the receivers.

Slide the padding for the grabrail onto the rail.

### 5.9 Hook and loop foot straps with instep pads

Attach the foot straps to the heel part of the footplate with two rivets or screws positioned beside one another (Fig. 39). In addition, the foot straps may also be used for the forefoot area.

Simply open and close the foot straps with the hook and loop closure.

### 5.10 (⊗) Headrest

Two headrest models are available; one for the use with and one for the use without padding.

When using the headrest without padding (Fig. 40) run the headrest straps around the rear of the backrest. Connect the buckle and pull the straps tight.

The headrest, with matching padding (Fig. 43), is screwed to the center part of the back padding at the desired height (Fig. 47). Both are removed together before the buggy/ stroller is folded.

### 5.11 Swivel lock

The swivel lock is installed before delivery by your authorized dealer or the manufacturer.

To open the swivel lock, turn the locking bolt by 90° until the bolt engages in the bore hole of the caster fork.

Disengage the swivel lock by pulling the locking pin out of the caster fork and turn the locking pin by 90° so that it remains in open position.

### 5.12 (⊗) *Padding*

The middle part of the padding (Fig. 43) is provided with an insert, to which the headrest, shoulder straps of the five-point belt and groin strap are connected. Therefore, it is possible to remove "everything with one grip", before folding the buggy/ stroller.

To mount the padding, you proceed as follows.

- Place the padding onto the seat.
- Pull the holding strap (Fig. 44) underneath the backrest, wrap the strap around the seat edge and affix the hook and loop straps.
- Wrap the flap at the rear side of the padding (Fig. 45) around the upper edge of the backrest and affix it to the hook and loop closure.
- Pull the side paddings (Fig. 46) around the left and right sides of the front frame (the armrests, if possible) towards outside and affix the hook and loop closure.
- To remove the padding, please proceed in the reverse order.

The back padding has zippers with double zipper pulls to provide openings (Fig. 47) for the headrest or the shoulder straps of the five-point belt.

### 5.13 (⊗) *Rain cape (Fig. 48)*

With the zipper open pull the hood over the child's head, then unfold the upper brim over the backrest and the lower brim around the footrest.

### 5.14 (⊗) *Slipsacks*

Place the summer, winter or sheepskin slipsack in the buggy/ stroller (Fig. 49). Tie the upper straps to the push handles. Tie the lower straps (Fig. 50) to the lower folding levers of the back. If using an abductor or belts, make a cut in the corresponding synthetic leather strips using a sharp blade.

### 5.15 (⊗) *Lateral guides with variable seat width reduction*

The pads of the lateral guides (Fig. 51) are filled with foam inserts of 1 cm ( $\frac{1}{3}$ " ) thickness which you can use variably. Thus the seat width can be reduced by 4, 6 or 8 cm (1.5", 2", 3").


To insert the lateral guides, you first plug the lateral guide padding (Fig. 52) with the insert onto the connection between the backrest and the frame and then you press down the insert onto the connection between seat tube and front frame.

### 5.16 (⊗) *Groin strap*

To fit the groin strap, you proceed as follows.

- Position the groin strap (Fig. 53) with the quick-release buckles on the thigh belts pointing downwards on the seat surface
- Thread the attachment tabs of the groin strap between seat and back padding (Fig. 54), wrap them around the seat edge and fix them using the hook and loop fasteners. If using padding, affix to the padding with the hook and loop fasteners.

- Thread the tensioning straps for the thigh belts on the left and right side through the eyelet (Fig. 55) and the loop end of the tensioning belt. Then tighten the tensioning strap.
- Thread the straps into the buckles (Fig. 56).

- You can now seat the child into the buggy/stroller on the groin strap. Guide the thigh belts (Fig. 57) up between the legs and over the thighs to the outside. The part of the thigh belt to which the buckle is attached should lie on the pelvic bone.
- Insert the male side of the quick-release buckles (Fig. 58) located on the tensioning straps into the buckles on the thigh belts.
- Tighten the thigh belt by pulling the free end of the strap that comes out of the buckle.
- To release the quick-release buckle, press the tab on the plug.

## 6 Using your LISA for transportation in a motor vehicle.

### **Danger!**

#### **Using your Ottobock Product as a Seat for Transportation in Motor Vehicles for Disabled Passengers**

We recommend that, wherever and whenever possible, users of a wheelchair, mobility base with seating shell or buggy transfer to the seats installed in the motor vehicle and use the corresponding vehicle restraint systems, because this is the only way to ensure optimal protection of the passengers in case of an accident.

**Your Ottobock product Lisa is permissible for transporting passengers in motor vehicles when using Ottobock safety components as well as appropriate restraint systems.** For more information, please refer to our instructions for use manual “**Using your wheelchair/mobility base with seating shell or buggy/stroller as a seat for transportation in motor vehicles**”, order no. 646D158.

## 7 Maintenance and Cleaning

### **Danger!**

The correct function of the rehab folding buggy/ stroller should be checked before each use.

Your LISA Rehab Folding Buggy/Stroller has been provided with the CE-marking. The manufacturer herewith guarantees that this medical product as a whole meets the requirements of European Directive 93/42/EEC.

The correct function of the LISA should be checked before every use.

The items listed in the following table must be checked by the user at the indicated intervals. Failure to carry out these simple checks may lead to problems arising that could invalidate the warranty.

| Check  | Before every use | Weekly | Monthly |
|--|------------------|--------|---------|
| Function test of the brake/wheel lock                                | X                |        |         |
| Check of the screw connections                                       |                  | X      |         |
| Air pressure (printed on the sidewall of the tire)                   |                  | X      |         |
| Visual examination of wearing parts, for example wheels and bearings |                  |        | X       |
| Contamination on bearings  |                  |        | X       |

Should any defects become apparent, please contact your authorized dealer to eliminate them.

We also recommend that you have the device serviced every twelve months by your Ottobock approved dealer.

**Instructions for Care**

- Clean all frame components and plastic parts using mild detergents only.
- Padding parts can be washed at 40 °C. If washed in a washing machine, put them in a linen bag or a pillow case.
- In most cases, wiping with a damp cloth is sufficient.
- Do not use salt water.
- Keep sand or other particles from damaging the wheel bearings.

## 8 Technical Data

| Dimensions (cm) and weights (kg)                | Size 1                             | Size 2                             |
|---|------------------------------------|------------------------------------|
| Model Variants                                  | HR32130000                         | HR32140000                         |
| Push handles (height of the buggy/ stroller)    | 103 (41")                          | 103 (41")                          |
| Overall width                                   | 60.5 (24")                         | 65.5 (26")                         |
| Min./max. overall length                        | 102 / 122 (40"/48")                | 102 / 122 (40"/48")                |
| Front wheel diameter                            | 18 (7")                            | 18 (7")                            |
| Rear wheel diameter                             | 25 (10")                           | 25 (10")                           |
| Folding size (LxWxH)                            | 103 x 45 x 43<br>(41" x 18" x 17") | 103 x 45 x 43<br>(41" x 18" x 17") |
| Turning radius                                  | 140 (55")                          | 145 (57")                          |
| Seat depth                                      | 30 / 35 (11.79"/13.79")            | 35 / 40 (13.79"/15.75")            |
| Seat width (without/with padding)               | 36 / 32 (14 1/5"/13")              | 41 / 37 (16/14.5")                 |
| Lower leg length                                | 27 - 40 (10.5"-15.75")             | 27 - 40 (10.5"-15.75")             |
| Seat height anterior/posterior                  | 50 / 38 (19 7/10"/15")             | 52 / 38 (20.9"/15")                |
| Backrest height                                 | 65 (25 3/5")                       | 70 (27,5")                         |
| Inclination of backrest towards the seat bottom | 90° to 120°                        | 90° to 120°                        |
| Seat angle                                      | 22° (110 lbs)                      | 22° (132 lbs)                      |
| Load capacity                                   | 50 (30 lbs)                        | 60 (31 lbs)                        |
| Weight  | 13.5                               | 14                                 |

**Note!**

Accessories and add-on components will reduce the remaining load capacity for the user.

## 9 Re-use / Disposal

### 9.1 Information on Re-use

**Danger!**

Risk due to improper use. Seat paddings which come into direct contact with the skin may cause functional or hygienic risks for another patient using one and the same product. They should be replaced in case of re-use.

The product is suitable for re-use.

Similar to second-hand machines or cars, products that are being re-used are subject to increased strain. Features and functions must not change in a way that could endanger patients or other persons within the product's life cycle. Based on market observations and the current state of technology, the manufacturer has calculated that the product drive can be used for a period of **4 years**, provided that it is used properly and that the service and maintenance instructions are observed. Periods during which the wheelchair is stored at the dealer or with the cost bearer are not included in this period. It should be clearly pointed out, however, that the product is reliable far beyond this defined period of time, provided that it is cared for and maintained appropriately. In cases of re-use, the corresponding product must first be thoroughly cleaned and disinfected. Afterwards, the condition of the product must be examined by an authorized dealer to check for wear and tear as well as any damage. Any worn and damaged components as well as components which do not fit or are unsuitable for the new user must be replaced.

The service manual includes a service schedule, detailed information, and a list of the required tools.

### 9.2 Information on Disposal

If a product is to be disposed of, all components of the product must be disposed of properly in accordance with the applicable environmental protection regulations in the respective country.

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Otto Bock Korea HealthCare Inc.
Otto Bock South East Asia Co., Ltd.
Otto Bock HealthCare GmbH