
ELLA STANDER



SN

TABLE OF CONTENT




1. Intended Use
2. Conformity
3. Terms of Warranty
4. Important Information
5. Indications, Contraindications and Application Risks
6. Stander Set Up (Assembly)
7. Initial Setting for User
8. Adjustments in Daily Use
9. Cleaning and Care
10. Product Inspections
11. Re-issuing of the ELLA Stander
12. Product Servicing
13. Technical Information
14. Product records
15. Distributor

Package content

- ELLA Stander according to ordered configuration.
- Charger (CH01) for battery when powered version is ordered.

Symbols used in Labels

This product label is an illustrative example used in explaining the symbols.



Product identifier
 9310: Size 1, Manual Tilt
 9320: Size 2, Manual Tilt
 9320+: Size 2+, Manual Tilt
 9330: Size 3, Manual Tilt
 9311: Size 1, Powered Tilt
 9321: Size 2, Powered Tilt
 9321+: Size 2+, Powered Tilt
 9331: Size 3, Powered Tilt
 9341: Size 4, Powered Tilt



CE-mark symbol



Medical device symbol



Date of manufacture
Year-Month-Day



Serial number




Maximum user weight
Kilograms



Manufacturer
Name and address



Consult instructions for use



Instructions for use translated by...
Translator contact information



Unique device identification number

Power Unit Labels



1. INTENDED USE

ELLA stander are mainly used with children, teenagers and young adults suffering of for example Cerebral Palsy (CP), muscular dystrophy, spina bifida and/or spinal cord injuries. It is possible to use ELLA Standers in any similar cases where the rehabilitation of individual requires or individual benefits of assisted standing. The conditions where ELLA Stander is used are to be determined by specialized health care experts such as doctors and physiotherapists.

ELLA Stander is a module based stander. Each module comes in different sizing to accommodate perfectly to end user body measures and creating a stable posture. All modules in any size can be combined with any ELLA stander size. ELLA Standers are available in four sizes. Sizes available can be found from "Technical Information" sheet under section 13. ELLA Standers are available with manual and powered tilt. The essential performance of the tilt is to safely transfer the end user between horizontal and vertical positions.

With ELLA standers it is possible to perform standing therapy in upright, prone and supine postures. All postures are possible to apply with or without hip abduction. Switching between postures is made easy and no additional modules or components is needed when switching from supine to prone and back. ELLA Standers are designed so that all its functions can be operated by using one 5mm hex tool or fingers. Tool is included to Stander.

ELLA Standers are designed in a way that all adjustments of different functions can be made outside the stander by using the 5mm hex tool or fingers. There is no need to remove any parts or cushions off the stander while doing adjustments.

All ELLA Stander functions are bilateral, individually adjustable and stepless. This allows the end user standing posture to be set independently for both right and left side. This feature enables the optimal set up for example to different leg lengths.

2. CONFORMITY

Declaration of Conformity

Timoteos Oy as manufacturer declares that the ELLA Stander as Class I medical device conforms to the requirements of EU Directive 93/42/EEC and the harmonized standard EN 12182:2012 when used as intended by the manufacturer.

Electrical Safety and EMC

This product complies with EN 60601-1:2007 General Requirements for basic safety and essential performance and EN 60601-1-2:2015 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests.

3. TERMS OF WARRANTY

A two year warranty applies to Timoteos products and modules from the date of purchase from Timoteos Oy. The warranty applies only when the product is used according to the specified conditions and for the intended purposes, following manufacturer's recommendations.

Availability of Spare Parts

Timoteos will supply spare parts to its products for 5 years after the product is manufactured. Manufacturing date is shown on the product label.

4. IMPORTANT INFORMATION

ELLA Stander is a multi-adjustable and versatile standing device for children, teenagers and young adults, that helps to achieve a specific and individual, comfortable standing posture. Product is to be used only in situations described in chapter 1 "Intended Use".

- 4.1. Always read instructions fully before use.
- 4.2. Ensure the product is fully assembled following instructions before use. Contact your local distributor if need advice at this stage.
- 4.3. Users should not be left unattended at any time whilst using Timoteos product.
- 4.4. Do not allow other individuals to use, play with, stand on, ride on or run around the Stander.
- 4.5. Only use Timoteos approved components with your product. Never modify the product in any way. In case modification is needed always consult Timoteos before making modifications. Service or repair should only be carried out by a Timoteos authorized service person. All services, inspections and possible modifications should be recorded. Failure to follow instructions may put the user or carer at risk and will invalidate the warranty of the product. Timoteos cannot accept any liability for injury or damage caused by unauthorized modification or unapproved components.
- 4.6. If in any doubt as to the continued safe use of your product or if any parts should fail, please cease using the product and contact Timoteos directly or your local distributor immediately.
- 4.7. Carry out all positional adjustments and ensure that they are securely fastened before you put the user into this product. Some adjustments may require the use of a tool which is provided with each product. Keep all tools out of reach of children and other people who has no need for using the tool.

- 4.8. All positional adjustments should be carried out by therapist or other technically competent person who has been trained to use of this product and has received instructions for the correct settings for individual user from the therapist. Timoteos recommends that a written record is maintained of all people trained to the use of this product.
- 4.9. This product could have been prescribed by the therapist or doctor who should have adjusted it to meet the user needs following the prescription. In such cases, if the product has not been adjusted by such person, please do not use the product and contact to therapist or doctor in question.
- 4.10. When placing the user into the ELLA Stander, for safety reasons, always secure the pelvic and trunk belts first. The knee supports should then be fastened and lastly secure the foot straps. When used in any position it is important to ensure that all belts and straps are fastened securely before making other adjustments or angling the stander to upright position. Always check that Buckle clips and Velcro's are fully engaged. If the user is wearing sandals or light footwear, short sleeve shirts or short sleeve trousers check the straps to make sure the webbing does not irritate the skin.
- 4.11. Always check with your therapist as to the correct tension of the straps for a particular end user as too much force could cause injury. Always use all positioning straps provided.
- 4.12. Fine positional adjustments may be carried out safely when the user is in the product. It is important to support all pads, straps and belts when adjusting them while user is in the Stander.
- 4.13. When doing a hoist transfer, please ensure that the hoist sling, or any part of hoist, is not caught on any part of the stander.
- 4.14. Always keep the tilt mechanism secured with tilt locking pin when tilt angle is not changed.
- 4.15. Although the ELLA Stander is fitted with castors it is not a mobility device. Always ensure that the castor brakes are locked at all times when the Stander is in use, being adjusted or even just stored.
- 4.16. When using the ELLA Stander, especially when tilting the Stander, ensure that the user and all parts of the product are well clear of surroundings furnishings to avoid potential collisions.
- 4.17. Never leave the product on a sloping surface. Always remember to lock all the castors.

- 4.18. The product contains components that are removable even not intended to remove and which could present a choking hazard to small children. Always check that knobs, bolts, washers and wingnuts are in place, tightened and secured at all times.
- 4.19. Product should be kept away from all direct sources of heat including naked flames, cigarettes, open fires, electric and gas heaters etc.
- 4.20. Never place hot items on the Tray as they may damage the plastic and can put the user or carer at risk of burn.
- 4.21. Never use the Tray to steer or push the Stander.
- 4.22. Ensure the Tray is always in place when using the Stander in Prone. Tray is for the use of the user only. Do not lean or place heavy objects on the Tray.
- 4.23. Clean the product regularly. Do not use abrasive cleaners. Carry out maintenance checks on a regular basis to ensure your product is in good working condition.
- 4.24. Do not carry out any maintenance or service to the product while product is in use.
- 4.25. Product is intended for indoor use in a home healthcare environment. When not in use product should be stored in dry place.
- 4.26. Manufacturer Product Serial Number can be found from the mobile base. Refer to that number in any correspondence. Product Serial Number is used to identify a single product and to track the content of the product configuration.
- 4.27. Before using ELLA Stander with electrical actuator for the first time, please charge the battery for 12 hours. Batteries should be charged overnight once every week to maintain optimum life of battery.
- 4.28. Do not operate the electrical tilt while charging the battery to avoid potential EM disturbance.
- 4.29. **WARNING!** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally. Use of accessories, transducers and cables other than those specified or provided by Timoteos could result in increased EM emission or decreased EM immunity of this equipment and result in improper operation. A power plug is supplied with this equipment. Please use the power plug suitable for your country. In case the supplied power plug is not suitable to your country, contact your local distributor for advice. Portable radio frequency (RF) communication equipment

(including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ELLA Stander, including cables specified by Timoteos. Otherwise, degradation of the performance of this equipment could result. This equipment is not likely to cause any harmful interference in nearby equipment; however, it can be impaired by electromagnetic (EM) interference emitted by other equipment, which can lead to malfunctions. Do not use this equipment where there is a risk that other electrical and/or electronic equipment can interfere with this equipment.

4.30. The safety advice and directions in this manual cannot cover all eventualities and unforeseen situations that may occur. It must be understood that common sense, general caution and care are factors that cannot be built into this product. These are factors which the person using this product must assume responsibility for. In case the use of this product causes a severe incident this should be immediately informed to manufacturer of the product in writing as well as to the responsible institution in your country following the guidelines given by this institution.

4.31. When a product reaches the end of its lifetime, it should be separated by material type so that the parts can be recycled or disposed of properly. If required, please contact your local dealer for a precise description of each material. Please contact your local authorities to be informed of the options available for recycling in an environmental friendly way.

5. INDICATIONS, CONTRAINDICATIONS AND APPLICATION RISKS

Indications

ELLA standers are used with children, teenagers and young adults suffering Paralysis of the lower extremities or muscle diseases, Cerebral Palsy (CP), muscular dystrophy, spina bifida and/or spinal cord injuries. It is possible to use ELLA Standers in any similar cases where the rehabilitation of individual requires or individual benefits of assisted standing. The conditions where ELLA Stander is used are to be determined by specialized health care experts such as doctors and physiotherapists.

Contraindications

Never use this product in any other condition than prescribed and determined by specialist healthcare expert.

In case user knee contracture is so severe that posture cannot be classified as standing posture, the product should not be used. The severity of each individual user knee contracture is to be determined by specialist healthcare expert.

To be able to use the product user must have adequate upper body control, such as head control, when product is used in Prone. Whether user has adequate upper body control is to be determined by specialist healthcare expert. In case user does not have adequate upper body control the product should be used in Supine posture only.

Application Risks

Never leave the user unattended to this product and do not allow other individuals to use, play with, stand on, ride on or run around the product. Failing to follow this may create a risk of injury to user or to other individuals.




This product is not intended to be a mobility device. Do not move the product while user is using the product. Always make sure that all the castors are locked when product is in use, adjusted or stored. Failing to follow this may create a risk of injury to user or to other individuals.

This product is to be used only indoors on an even surface. Failing to follow this may create a risk of product falling over.

6. STANDER SET UP (ASSEMBLY)

STANDER SET UP

3-in-1 and 4-in-1 Stenders



SETTING UP THE STANDER

- Remove the protective filling from the box and take the base (1) and frame (3) out of the box. Base is folded in the box, in Frame all the ordered modules are already attached.
- Lift the tilt / Frame holder up (2), place the two bolts into their slots and tighten the bolts to secure the tilt / Frame holder to its position. If the tilt / Frame holder does not slide into its place properly, loosen the bolts as long as tilt/ Frame holder slides into its place. It is not needed to remove the bolts off the base.
- Place the Frame on the Base (3) by placing the three bolts into their slots and tighten the upper bolt to secure the Frame locking on the Base. If the Frame does not slide into its place properly, loosen the upper bolt as long as Frame slides into its place. It is not needed to remove the bolts off the Frame.
- ELLA Stander is now ready (4) for initial, user specific, settings.

7. INITIAL SETTING FOR USER



INITIAL SETTING FOR USER

60-72 cm in ELLA Size 1 and 90-105 cm in ELLA Size 2
3-in-1 and 4-in-1 Standers

NOTICE

Following instruction, on this page apply only to ELLA 1 in user height range of 60-72 cm and ELLA 2 in user height range of 90-105 cm.

Because of design that allows wider growth range the starting point of the initial setting of ELLA 1 and ELLA 2 for small users is determined by the lowest tray table height where as normally the starting point is lowest foot plate height.



SETTING UP THE STANDER

- Place the tray holder and tray to its place in the lowest possible position (1).
 - Place the lower edge of the trunk lateral support (in its lowest position) to the same height with tray table top by adjusting the height of the trunk support.
 - Measure the width of the users trunk (2) and adjust the trunk lateral support to equal width.
 - Measure the height of the user from sole to armpit (3). Adjust the foot plate up or down so that measure from bottom of the foot plate to top edge of the trunk lateral support is equal.
 - Measure the height of user from sole to middle of pelvis (4). Adjust the pelvic support pad so that the measure from the bottom of the foot plate to middle of the pelvic support pad is equal.
 - Measure the width of the users hip (5) and adjust the pelvic lateral supports to equal width.
 - Measure the height of user from sole to right below patella (6). Adjust the knee support up or down so that the measure from the bottom of the sole to the upper edge of the knee support is equal.
 - Adjust Head Support to right height and angle for user.
- Notice that in ELLA Stander the user standing posture is possible to finalize when user is placed on stander and stander tilt in is in desired angle.


INITIAL SETTING FOR USER

All user heights in all 4 stander sizes except the ones instructed on previous page
3-in-1 and 4-in-1 Stander

NOTICE

Following instruction, on this page, apply to ELLA 1 in user height range of 73-95 cm, ELLA 2 in user height range of 106-130 cm and all user heights in ELLA 3 and ELLA 4.

Because of design that allows wider growth range the starting point of the initial setting of ELLA 1 and ELLA 2 for small users is determined by the lowest tray table height where as normally the starting point is lowest foot plate height. Set up for small users was presented on earlier page.



SETTING UP THE STANDER

- Place the foot plates (1) in to lowest position.
 - Measure the height of user from sole to right below patella (2) or right above patella in case knee support is wished to put on top of patella. Adjust the knee support up or down so that the measure from the bottom of the sole to the upper edge (or lower edge in case support will come above patella) of the knee support is equal.
 - Measure the height of user from sole to middle of pelvis (3). Adjust the pelvic support pad so that the measure from the bottom of the foot plate to middle of the pelvic support pad is equal.
 - Measure the width of the users hip (4) and adjust the pelvic lateral supports to equal width.
 - Measure the height of the user from sole to armpit (5). Adjust the trunk support pad up so that measure from bottom of the foot plate to top edge of the trunk lateral support is equal.
 - Measure the width of the users trunk (6) and adjust the trunk lateral support to equal width.
 - Adjust Head Support to right height and angle for user.
- Notice that in ELLA Stander the user standing posture is possible to finalize when user is placed on stander and stander tilt in is in desired angle.

4-in-1 ADDITIONAL SETTINGS

Only 4-in-1 Standers

NOTICE

The following instructions apply only to ELLA standers which are bought as 4-in-1 or which has been upgraded with 4-in-1 Assembly Kit.

1)



Setting the angle for abduction

- Download a free clinometer app or use manual clinometer.



- Pre-Set the angle when stander frame is in vertical position by placing a bar or stick to simulate leg and measuring the angle with clinometer (2).

➤ Notice that it is possible to adjust both legs individually even to different angles



- Once user is standing and set up is finalized, double check the angle by measuring it with clinometer from the user leg. Correct angle if needed.

➤ Notice that small finalizing changes in angle are possible to do while user is in stander and stander frame is in vertical position. If the foot module is difficult to move, adjust the stander frame towards horizontal position to ease the pressure from the foot allowing the foot module to slide more easily.

PUTTING THE USER ON STANDER

3-in-1 and 4-in-1 Standers

NOTICE

It is recommended to use lift or hoist when placing the user to stander. If this is not possible it is advised to pay attention to user weight and proper ways of lifting to avoid injuries.

It is recommended to put the user on stander while stander is in horizontal position.

Always make sure when user is in the stander that all the belts and supports are in place and tightened. Do not leave child unattended to stander.



PUTTING THE USER ON STANDER

- Place the stander frame into horizontal position (1) and open all the straps and belts. It is advised to turn the belts and straps aside before placing the user on stander. Open and turn aside the possible swing-away lateral supports from the side where user will be lifted to stander.
- Notice there is small pieces of velcro to help the lateral support belts to stay aside.
- Fasten all the belts and straps and close all the swing-away supports
- Pay attention that belts on pelvis and trunk area are firmly in place and check there is no empty space between belt and body (2). This way the user posture in stander remains as stable as possible.
- Check that the knee support is correctly placed right under patella (3) or right above patella in cases where it is wished to put so.
- Lift the stander frame into desired angle (4).
- Notice it is normal that user slides downwards when angle of stander frame changes from horizontal towards vertical. This may impact to pre-set adjustments. Once the angle is set it is possible to finalize the adjustments of the modules while user is in stander and stander frame is in desired angle. Simply loosen the bolts of modules and finalize the location of the module.

Fastening the belts

3-in-1 and 4-in-1 Standers



BELT / STRAP TYPES

- 1) Buckled strap in Flexible Swing-Away Lateral Supports.
- 2) Velcro enhanced buckle straps with cover in Standard Lateral Supports and Swing-Away Lateral Supports.
- 3) Buckled strap in knee support
- 4) Velcro strap in knee support
- 5) Buckled strap with cover in Foot Rest (Sandals)



Flexible Swing-Away

- 1) Pre-adjust the length of straps.
- 2) Close the Flexible Swing-Away Laterals
- 3) Fasten the buckle and tighten if needed by pulling the strap outwards.



Knee Support

- 1) Pre-adjust the length of straps.
- 2) Turn the knee pads over the leg
- 3) Slide the velcro through metal ring and tighten or fasten the buckle and tighten if needed by pulling the strap outwards.



Trunk and Pelvic Belts

- 1) Pre-adjust the length of straps. Notice the small piece of velcro outside the belt.
- 2) Turn the belt outside and use the small piece of velcro to fasten it to lateral support. This helps to keep the belts away when placing the user to stander.
- 3) Wrap the belts around user, fasten the buckle and tighten if needed by pulling the straps outwards.
- 4) Check that there is no empty space between the belt and body. Otherwise the upper body posture might not be as stable.


8. ADJUSTMENTS IN DAILY USE

Adjustment locations

- A. Knee Support height, depth, width and rotation adjustment.
- B. Lateral Support (trunk and pelvic) height, width and rotation adjustment.
- C. Head Rest rotation adjustment.
- D. Head Rest depth adjustment.
- E. Head Rest height adjustment.
- F. Trunk Support Pad height adjustment.
- G. Tray height adjustment.
- H. Pelvic Support Pad height adjustment.
- I. Gas Piston assisted tilt adjustment pedal.
- J. Gas Piston Assisted tilt locking pin. Always keep the tilt locked when angle is not changed.
- K. Castor locking pedal.
- L. Foot Rest height adjustment.
- M. Foot Rest width adjustment.
- N. Foot Rest depth adjustment.
- O. Rotational Trunk and Pelvic Support adjustment (2 bolts on both sides)




- P. Handset HB-41 for electrical tilt operating from horizontal to vertical
- Q. Charging port for battery charger.
Do not operate the electrical tilt while charging the battery.
- R. Electrical tilt emergency stop




Swing-Away Laterals (including Flexible version)

1. Pull the knob up to release lock of Swing-Away Lateral Supports
2. Move the Lateral Support sideways.




Shoulder Protraction Pads

1. Loosen the wingnut to release the lock of Protraction Pad. Move Pad to desired location and tighten the wingnut.



Harness Guides

1. Loosen the wingnut to release the lock of Harness Guide. Move the Guide sideways to desired width and tighten the wingnut. Adjust both sides separately.



Change from Supine to Prone

1. Remove the Head Rest by loosening head rest height adjustment knob and pulling the head rest off the stander.
2. Insert the tray from the side of tilt mechanism (opposite to supine set up).
3. Turn the Sandals 180° by loosening the bolt underneath the Sandals.



9. CLEANING AND CARE

When cleaning the product, we recommend that you use only warm water and a nonabrasive detergent. Never use organic solvents or dry-cleaning fluids.

Staining should be removed from upholsteries as quickly as possible with absorbent cloth, towels or a sponge. Soap and warm water sponging are effective for ordinary soiling and minor spills. Do not over wet the fabric, as this will cause the staining to spread.

Always ensure the product is dry before use.

Always replace broken or damaged upholsteries and cushions. Upholsteries and cushions are available as spare parts and we recommend contacting your local distributor for new ones.

For deep cleaning of upholstery, a low-pressure steam cleaner can be used in situ. Machine washable upholstery is available as optional accessory. Contact your local distributor for machine washable upholstery.

Soap and water or antibacterial spray can be used for daily cleaning of metal and plastic components. For deep cleaning a low-pressure steam cleaner can be used.

Do not use solvents to clean plastic or metal components.

Always ensure product is dry before use.

Timoteos recommends the product is stored and transported in the temperature range of -10°C to +50°C, under relative humidity of 20% to 90% and under atmospheric pressure of 800 to 1060 hPa. Operating temperature range is +5°C to +40°C, relative humidity and atmospheric pressure as above.

Avoid storing the product for long time under extreme temperatures, temperature variations, in high humidity or in direct sunlight.

10. PRODUCT INSPECTIONS

We recommend that daily visual inspections of the product are carried out by therapist, carers or parents to ensure the product is safe for use. The recommended daily inspections are minimum as detailed below.

- Ensure all adjustment knobs and bolts are in place, tightened and secured.
- Check all upholstery for signs of wear and tear.
- Check all wheels are moving freely and lock securely.
- Check all Velcro straps and brush fluff to ensure secure firmly.
- Check all buckle straps to ensure the buckle functions and secure firmly.
- Check that all electrical connectors, components and cables are intact.

Timoteos recommend that each product should be subjected to a detailed inspection as least once a year and every time the product is reissued for use or any part is changed or replaced. This inspection should be carried out by a technically competent person

who has been trained in the use of the product and should include the following checks as a minimum requirement.

- Check all knobs, nuts, bolts and plastic buckles are in place. Replace any missing items. Paying particular attention to the following items;
 - ✓ Buckles and Velcro bands and straps on pelvic and trunk, knee support straps, buckles and Velcro, and foot rest buckles and straps.
 - ✓ Trunk, Hip, Knee and Foot supports adjust and locks securely.
 - ✓ Head Rest adjust and locks securely.
- Lift the base and check each castor individually. Make sure they are moving freely and remove any dirt from the wheels. Check that the brakes lock the castors securely.
- Check that tilt mechanism functions and the gas piston works properly. Check that the tilt locking pin locks securely.
- Visually check the structure of the product ensuring there are no signs of fatigue or cracking around the welds.
- Visually check the plastic parts of the product ensuring there are no signs of fatigue or cracking.
- Check that all electrical connectors, components and cables are intact.

Timoteos recommends that a written record is maintained of all product inspections.

11. RE-ISSUING OF THE ELLA STANDER

Most Timoteos products are assessed and ordered to meet the needs of an individual user. Before re-issuing a product we recommend that the therapist prescribing (or other medical expert entitled to prescribe) the product has carried out an equipment compatibility check for the new user and has ensured that the product being re-issued contains no modifications or special attachments.

Additionally a detailed technical inspection should be carried on the product prior to re-issuing. This should be carried out by a technically competent person who has been trained in the use and inspection of the product. Please refer to section 10 for the required checks to be carried out.

Before re-issuing ensure the product has been cleaned thoroughly in accordance with section 9.

Ensure a copy of the user instructions is supplied with the product. A copy is available from Timoteos upon request.

Timoteos recommends that a written record is maintained for product re-issuing.

12. PRODUCT SERVICING

Servicing of Timoteos products should only be carried out by technically competent persons who have been trained in the use of the product. All service enquiries should be directed to local Timoteos distributor who will assist you further. In case a local distributor is not available, please contact to Timoteos Oy directly for assistance. Timoteos recommends that a written record is maintained of all product services.

13. TECHNICAL INFORMATION

ELLA STANDING FRAME	TECHNICAL INFORMATION							
	Size 1		Size 2		Size 3		Size 4	
	min	max	min	max	min	max	min	max
User approximate age (years)	0	4	3	7	6	12	11	18+
User height [mm]	600	950	900	1300	1200	1550	1300	1800
Max user weight [kg]	30		40		55		80	
Power assisted tilt	Yes		Yes		Yes		Yes	
Manual assisted tilt	Yes		Yes		Yes		No	
Pivot width [mm]	1000		1000		1350		1450	
Castor height (space for hoist skis) [mm]	140		140		140		140	
Castor diameter [mm]	100		100		100		100	
Height in upright position, floor to top of head rest [mm]	970	1100	1040	1430	1280	1630	1460	1880
Height in horizontal position to lateral support edge [mm]	810		890		940		950	
Height in horizontal position to trunk support pad surface [mm]	710		790		790		810	
Foot plate to top of head rest [mm]	600	850	740	1300	1040	1600	1075	1800
Foot plate to top of trunk pad [mm]	430	680	570	975	800	1280	900	1500
Foot plate to top of pelvic pad [mm]	300	540	475	785	500	1040	625	1225
Foot plate to top of knee cup [mm]	140	450	180	550	200	650	250	850
Foot plate to floor [mm]	135	-	135	-	150	-	150	-
Product weight with manual tilt, Standard Set up [kg]	28		29		32		-	
Product weight with powered tilt, Standard Set up [kg]	32		35		38		47	
Width [mm]	580		580		702		752	
Length [mm]	750		750		1000		1080	
Actuator Voltage (LA31, LA34)	24V (LA31)						24V (LA34)	
Integrated Lead Battery / Control Unit (CB8)	24V / 2,9 Ah (Control Box CB8)							
Operating noise with powered tilt	48 dB (LA31)						< 50 dB (LA34)	
Lifting Speed without / with load [mm/s]	8,2 / 5,4 (LA31)						10,8 / 7,9 (LA34)	
Battery charger (CH01)	Mains: 100 - 240 VAC / 50 - 60 Hz switch mode supply Charging Voltage: 27,6 VDC Charging Current: Max. 500 mA							
Battery charging time	6 hours							
IP Protection Class	IPX1 (Handset HB41) IPX4 (Actuator LA31, LA34)							
Duty Cycle	10%, Max. 2 min / 18 min (LA31, LA34, CB8)							
Emergency stop	Yes							
Power tilt components	ELLA 1 -3: Actuator LA31, Control Box CB8, Charger CH01, Handset HB41 ELLA 4: Actuator LA34, Control Box CB8, Charger CH01, Handset HB41							

ELLA STANDING FRAME	TECHNICAL INFORMATION							
	Size 1		Size 2		Size 3		Size 4	
	min	max	min	max	min	max	min	max
User approximate age (years)	0	4	3	7	6	12	11	18+
User height [in]	23,62	37,40	35,43	51,18	47,24	61,02	51,18	70,87
Max user weight [lb]	66,14		88,18		121,25		176,37	
Power assisted tilt	Yes		Yes		Yes		Yes	
Manual assisted tilt	Yes		Yes		Yes		No	
Pivot width [in]	39,37		39,37		53,15		57,09	
Castor height (space for hoist skis) [in]	5,51		5,51		5,51		5,51	
Castor diameter [in]	3,94		3,94		3,94		3,94	
Height in upright position, floor to top of head rest [in]	38,19	43,31	40,94	56,3	50,39	64,17	57,48	74,02
Height in horizontal position to lateral support edge [in]	31,89		35,04		37,01		37,4	
Height in horizontal position to trunk support pad surface [in]	27,95		31,1		31,1		31,89	
Foot plate to top of head rest [in]	23,62	33,46	29,13	51,18	40,94	62,99	42,32	70,87
Foot plate to top of trunk pad [in]	16,93	26,77	22,44	38,39	31,5	50,39	35,43	59,06
Foot plate to top of pelvic pad [in]	11,81	21,26	18,7	30,91	19,69	40,94	24,61	48,23
Foot plate to top of knee cup [in]	5,51	17,72	7,09	21,65	7,87	25,59	9,84	33,46
Foot plate to floor [in]	5,31	-	5,31	-	5,91	-	5,91	-
Product weight with manual tilt, Standard Set up [lb]	61,73		63,93		70,55		-	
Product weight with powered tilt, Standard Set up [lb]	70,55		77,16		83,78		103,62	
Width [in]	22,83		22,83		27,64		29,61	
Length [in]	29,53		29,53		39,37		42,52	
Actuator Voltage (LA31, LA34)	24V (LA31)						24V (LA34)	
Integrated Lead Battery / Control Unit (CB8)	24V / 2,9 Ah (Control Box CB8)							
Operating noise with powered tilt	48 dB (LA31)						< 50 dB (LA34)	
Lifting Speed without / with load [in/s]	0,32 / 0,21 (LA31)						0,43 / 0,31 (LA34)	
Battery charger (CH01)	Mains: 100 - 240 VAC / 50 - 60 Hz switch mode supply							
Battery charging time	6 hours							
IP Protection Class	IPX1 (Handset HB41) IPX4 (Actuator LA31, LA34)							
Duty Cycle	10%, Max. 2 min / 18 min (LA31, LA34, CB8)							
Emergency stop	Yes							
Power tilt components	ELLA 1 -3: Actuator LA31, Control Box CB8, Charger CH01, Handset HB41 ELLA 4: Actuator LA34, Control Box CB8, Charger CH01, Handset HB41							

Testing for EM disturbances

Test modes during emission tests

1. Battery mode: Standby and moving
2. Use when charging: Standby

EM emission standards	Classification / levels
Conducted and radiated RF emissions CISPR 11	Group 1 Class B
Harmonic distortion IEC 61000-3-2	Class A
Voltage fluctuations and flicker IEC 61000-3-3	$d_{max} \leq 4\%$

Test modes during immunity tests

1. Battery mode: Standby
2. Use when charging: Standby

EM immunity standards	IEC 60601-1-2:2014 test levels	
Electrostatic discharge IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	
Radiated RF EM fields IEC 61000-4-3	10 V/m 80 MHz - 2,7 GHz 80 % AM at 1 kHz	
Proximity fields from RF wireless communication equipment IEC 61000-4-3	Test frequency (MHz)	Test level (V/m)
	385	27
	450	28
	710	9
	745	
	780	
	810	
	870	28
	930	
	1720	
	1845	
	1970	
	2450	9
	5240	
5500		
5785		
Rated power frequency magnetic fields IEC 61000-4-8	30 A/m 50 Hz	
Electrical fast transients / bursts IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	
Surges Line-to-line IEC 61000-4-5	$\pm 0,5$ kV, ± 1 kV	
Surges Line-to-ground IEC 61000-4-5	$\pm 0,5$ kV, ± 1 kV, ± 2 kV	
Conducted disturbances induced by RF fields IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz	
Voltage dips IEC 61000-4-11	0% U_T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T ; 1 cycle 70 % U_T ; 25 cycles at 0°	
Voltage interruptions IEC 61000-4-11	0 % U_T ; 250 cycles at 0°	

14. PRODUCT RECORDS

Date	Action	Signature	Print name

15. DISTRIBUTOR