# Hurtigguide MEMOplanner

---

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Side 3 (9)

5. Vælg dato.
6. Vælg type.

The image displays three sequential screenshots from a mobile application for adding an activity.

**Left Screenshot: "Vælg dato"**  
Shows a calendar for October 2016. The date 27 (Tuesday) is selected. The calendar grid shows days from 1 to 31, with the week starting on Monday (U). The selected date is highlighted in pink.

**Middle Screenshot: "Ny aktivitet"**  
Shows a form for adding a new activity. The form has tabs for "Navn", "Ekstra", "Alarm", and "Påmindelse". The "Navn" tab is active. Fields include:

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Buttons at the bottom are "< Forrige" and "Næste >".

**Right Screenshot: "Vælg Type"**  
Shows a dialog for selecting the activity type. Options are:

- Hele dagen:
- Venstre:
- Højre:

Buttons at the bottom are "Annuller" and "OK".

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The image displays two screenshots from a mobile application. The left screenshot, titled 'Vælg tid' (Select time), shows a time selection interface with a numeric keypad. The start time is set to 16:00, and the end time is currently blank (--:--). The right screenshot, titled 'Ny aktivitet' (New activity), shows a form with several fields: 'Navn' (Name), 'Billede' (Image), 'Dato' (Date) set to 'tirsdag 27 september', 'Type' (Type) set to 'Højre', 'Tid' (Time), 'Kvitterbar' (Checkable), and 'Slet bagefter' (Checkable). Arrows point from the text on the right to the 'Kvitterbar' and 'Slet bagefter' checkboxes.

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