- Ⓛ **D** **Rollbarer Bett-Tisch**
- Ⓛ **GB** **Overbed table**
- Ⓛ **DK** **Sengebord**
- Ⓛ **I** **Tavolino da letto**

Art.Nr.
Item no. **272010**
Varenr.
Codice

MOBILEX A/S
Grønlandsvej 5
DK-8660 Skanderborg
Danimarca
www.mobilex.dk

D Rollbarer Bett-Tisch - Art. 272010

GB Overbed table - item no. 272010

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Rollable bed-side table with 4 wheels. The tray of the table can be easily adjusted in height and angle and is suitable to carry smaller items like books, food, laptop, tablet or similar). The weight of these items must not exceed the carrying capacity of the tray (see specifications below).

Safety requirements

- Make sure that this user manual is read by all persons using the device. Manufacturer doesn't take any responsibility for damages and/or injuries, caused by the fact that the user manual has not been followed.
- Only use the product in flawless condition.
- If defects or errors are detected you must immediately contact the dealer.
- Follow the information on the product label.
- Only use the product for the described purpose
- Do not exceed the carrying capacity of the tray

Product description

Rollable bed-side table with height and tilt adjustable wooden chipboard tray and a metal frame. It comes with 4 lockable wheels.

Assembly

Before assembly, you should check with reference to the drawing below that no parts are missing (parts A-D, 1-9). Then assemble the table following the number series on the drawing.

Cleaning & maintenance

Normal dirt on metal and plastic parts can be removed with standard cleaning agents. Check the specific product information and use only commercial cleaning products which are suitable for cleaning and disinfection (no solvents or abrasives).

Guarantee terms

Mobilex A/S offers 2 years of warranty for damages caused by errors in production or material errors. The warranty is not valid by not-intended repair or usage of the product. Parts subject to normal wear are not covered by the warranty, unless wear is caused by a manufacturing fault. If product changes are made without our written acceptance the CE-marking and product warranty will not be valid. In the event of product damages covered by this warranty, please notify your dealer or Mobilex A/S directly. The warranty does not cover transport costs and does not include compensation for personal injuries or non-use of the product during repair. The warranty does not cover damages due to the user manual not being followed.

Technical data

Overbed table	itemnr. 272010
Size of the tray (D x L)	41 x 55 cm
Height (adjustable)	from 72 to 115 cm
Weight	7,5 kg
Carrying capacity	15 kg

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