**VISION**  
by ropox

Vision

Tekniske data / Technical data / Technische Daten

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
Vision

Technical Data

**Keep always this folder with product!**

Contents:

<b>1.</b>	<b>INTRODUCTION</b> .....	<b>11</b>
<b>2.</b>	<b>COMPLIANCE WITH EU-DIRECTIVES</b> .....	<b>11</b>
<b>3.</b>	<b>APPLICATION</b> .....	<b>11</b>
<b>4.</b>	<b>ASSEMBLY</b> .....	<b>11</b>
<b>5.</b>	<b>USER GUIDE</b> .....	<b>11</b>
<b>6.</b>	<b>TECHNICAL DATA</b> .....	<b>12</b>
<b>7.</b>	<b>PERFORMANCE TEST</b> .....	<b>13</b>
<b>8.</b>	<b>SAFETY IN USE</b> .....	<b>13</b>
<b>9.</b>	<b>CLEANING/ MAINTENANCE</b> .....	<b>14</b>
	9.1 <i>Frame</i> .....	14
	9.2 <i>Tabletop, laminated surfaces</i> .....	14
	9.3 <i>Maintenance</i> .....	14
<b>10.</b>	<b>TROUBLE SHOOTING</b> .....	<b>14</b>
	10.1 <i>The table seems loose or unstable</i> .....	14
	10.2 <i>The table cannot be adjusted in height or tilt adjustment seems very difficult</i> .....	14
<b>11.</b>	<b>CE-MARKING</b> .....	<b>15</b>
<b>12.</b>	<b>GUARANTEE</b> .....	<b>15</b>

## 1. Introduction

VISION is the new range of modern height-adjustable Ropox tables. In our opinion a product must be functional and beautiful – no matter who is going to use it. That's what we call VISION.

Vision meets all needs and ensures ergonomically correct working positions for people with or without disability.

Vision is available in several sizes to suit individuals and carers needs. Sufficient legspace makes the table ideal for wheelchair users and gives the therapist optimum working conditions in a training situation. Due to the unique design with rounded corners there is no danger of hitting sharp edges when moving round the table and two inconspicuous castors at the back, make it easy to draw the table closer to the user. As an option extra arm supports may be clipped on and off as required.

To suit specific activities and individual requirements, the tables are easily and quickly adjustable in height from 50-70 cm or 60-90 cm by means of a detachable handle.

All tables are available in six colours: white, light grey, dark grey, red, blue and green with matching leg covers to the red, blue and green tabletops. Natural aluminium is used on the other colours.

This document should always accompany the product and be read by and be available to the user. Correct use, operation and inspection are decisive factors for efficient and safe performance.

As this product is adjustable in height there is a risk of squeezing. The product must therefore always be operated by or under the guidance of an experienced adult, who has read and understood the importance of section 8 "Safety in use".

It is very important that there are no obstacles between the tilted tabletop and the frame to prevent any change of entrapment.

## 2. Compliance with EU-Directives

This product has CE-marking according to the current Directives and thus complies with the basic safety requirements. See separate CE-declaration, section 11.

## 3. Application

Vision is designed to obtain optimum working height for the user. It is an activity table not intended for lifting people or objects. Vision should be used indoor, at temperatures and humidity as described in section 6.

## 4. Assembly

See separate assembly instructions.

## 5. User guide

See separate user guide.

## 6. Technical data

<b>Product:</b>	<b>Vision</b> Type A, B, C, D, E, F, G 90x60 cm, 120x60 cm, 90x90 cm, 165x100 cm, 200x100 cm, 240x100 cm
<b>Item numbers:</b>	Manual 20-6XXXX
<b>Height range:</b>	50–70 cm (manual) incl. tabletop 60–90 cm (manual) incl. tabletop
<b>Frame material:</b>	Aluminium leg. 606045-T6 Steel St. 37 Various plastic components
<b>Surface treatment:</b>	Anodic coating Powder coating, various colours Chromating
<b>Tabletop:</b>	Etronit-M laminate
<b>Max. load on the table (static):</b>	100 kg evenly distributed over the entire tabletop
<b>Max. load during lifting (dynamic):</b>	50 kg
<b>Max. load on the tabletop tilted (static):</b>	30 kg
<b>Max. load during tilting the tabletop (dynamic):</b>	4 kg
<b>Power supply:</b>	230VAC / 2.5A / 50 Hz - Alternatively 100VAC / 50/60 Hz
<b>Standby primary:</b>	1W
<b>Control voltage:</b>	24VDC
<b>Duty circle:</b>	Max. 10 % conform to 1 min. active / 9 min. pause
<b>Speed:</b>	Approx. 32 <sup>mm</sup> / <sub>sec</sub>
<b>Lifting time stroke 30 cm:</b>	Approx. 10 sec.
<b>Temperature:</b>	5-45°C
<b>Air humidity:</b>	5-85% (non-condensing)
<b>Guarantee:</b>	See guarantee, section 12
<b>Manufacturer:</b>	Ropox A/S, DK-4700 Naestved, Tlf.: +45 55 75 05 00 E-mail: info@ropox.dk, www.ropox.com

## 7. Performance test

After assembly and prior to use, all functions of the table must be tested. Subsequently this test must be carried out by competent personnel at least once a year:

1. Check that the assembly instructions have been observed.
2. Check that all bolts have been tightened.
3. Check that there is no load on the table.
4. Check that there is nothing preventing the table from moving freely within the height adjustment range.
5. Place the handle according to the “user guide” and move the table to the bottom position. Check that the movement is even and smooth.
6. Now move the table to the top position. Check that the movement is even and smooth.

### If Tilt fitting is mounted:

7. Place the handle for the tilt fitting according to the “user guide” and tilt the tabletop to the upright and down position. Check that the movement is even and smooth.

When the above mentioned items have been checked, the table is ready for use. See “Safety in use” section 8.

## 8. Safety in use

- The table should only be used by people who have read and understood these instructions.
- The table is not intended for lifting people or objects.
- Always use the table in a manner precluding the risk of damage to persons or property. The person, who operates the table, is responsible for avoiding damage or injury.
- Children and people with reduced observational ability must only be allowed to operate the table under supervision.
- When the table is used in publicly accessible locations where children or other people with reduced observational capability may get close to it, the person operating the table must pay utmost attention in order to prevent dangerous situations.
- Make sure that there is free space above and below the table to allow height adjustment.
- Do not operate the table in the case of errors or damage.
- Do not overload the table and make sure that the load distribution is correct.
- In connection with inspections, service or repairs make sure that the table is not loaded.
- Any modification of the table, which may influence its operation or construction, is forbidden.
- If the table has not been assembled in accordance with these mounting instructions, the guarantee may become void.
- Only use Ropox original spare parts as replacement parts. If other spare parts are used, the guarantee may become void.

## 9. Cleaning/ Maintenance

### 9.1 Frame

Clean the frame with a damp cloth using lukewarm water and ordinary cleaning agents. Do not use corrosive/ abrasive liquids or grinding cloths, brushes or sponges. Dry the frame after cleaning.

### 9.2 Tabletop, laminated surfaces

Clean the Etronit-M surface with a soft cloth firmly wrung in hot water with a little liquid detergent. The detergent should be without sulphonic soap as this may result in a greasy surface. Stains that cannot be removed with water can often be removed with a mixture of equal parts chlorine bleaching and water or with acetone. When using these products see that the room is sufficiently ventilated and note that these products should not be used on the frame.

Stains from felt-tip pencils can be removed using methylated spirit and pencil strokes can be removed with a piece of eraser. In case of doubt, we recommend you to experiment on a spot that cannot be seen as certain products may damage the laminate, e.g. sulphuric acid etc.

Some chemicals, ink and strongly colouring fruit juices may cause discoloration if they remain on the table for a longer period of time as the pigmentation sticks to the surface. Do not use steel wool, scouring powder and other abrasive products or detergents leaving a film on the surface, e.g. yellow soap.

After cleaning, wash the table with a soft cloth firmly wrung in clean water and dry with a cloth.

### 9.3 Maintenance

Inspections, service and repairs must be carried out by competent personnel.

The frame is maintenance-free and the moving parts have been lubricated for life. For reasons of safety and liability we recommend inspection of the frame once a year:

- Check that all bolts have been tightened.
- Check that the table moves freely and without problems from bottom to top position.

Only use Ropox original spare parts as replacement parts. If other spare parts are used, the guarantee may become void.

## 10. Trouble shooting

### 10.1 The table seems loose or unstable

- The screws assembling the frame have not been securely tightened.
- Tighten all screws of the frame. See assembly instructions.

### 10.2 The table cannot be adjusted in height or tilt adjustment seems very difficult

- Check that there is no load on the table
- Check the table for overload

## 11. CE-marking

### Vision, height-adjustable table

#### Tabletop 90x60 cm

Type	Height 50-70 cm	Height 60-90 cm
A	20-60010	20-60020
B	20-60011	20-60021
C	20-60012	20-60022
D	20-60013	20-60023
E	20-60014	20-60024
F		
G		

#### Tabletop 120x60 cm

Height 50-70 cm	Height 60-90 cm
20-60210	20-60220
20-60211	20-60221
20-60212	20-60222
20-60213	20-60223
20-60214	20-60224
20-60215	20-60225
20-60216	20-60226

#### Tabletop 90x90 cm

Type	Height 50-70 cm	Height 60-90 cm
A	20-60310	20-60320

#### Tabletop 165x100 cm

Height 50-70 cm	Height 60-90 cm
20-62116	20-62216

#### Tabletop 200x100 cm

Type	Height 50-70 cm	Height 60-90 cm
A	20-62120	20-62220

#### Tabletop 240x100 cm

Height 50-70 cm	Height 60-90 cm
20-62124	20-62224

conform to the following directives and standards:

### DIRECTIVES

EU-Directive 93/42EEC on Medical equipment, modified by Directive No. 98/79/EC

### STANDARDS

DS/EN	1041:	2000	DS/EN	12182:	2002
DS/EN	14971:	2001	DS/EN ISO	9999:	2002

Declarations of conformity can be supplied upon request.

## 12. Complaints


We refer to our general Terms of sale and delivery on our homepage [www.ropox.com](http://www.ropox.com)

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Konformitätserklärungen sind in Original bei Ropox A/S erhältlich.
## ROPOX A/S