



**PhysioNova**  
Health Care Systems



**GET OUT OF BED WITH THE ROTOFLEX<sup>®</sup>**  
**A STEP BACK INTO INDEPENDENCE**



**ROTOFLEX<sup>®</sup>**

Into and Out of Bed at the Touch of a Button.

# Get out of bed with the ROTOFLEX® – a step back into independence

---

**The familiar private surroundings of the patient are maintained.**

---

The **ROTOFLEX®** is an innovative nursing bed which was developed to improve the quality of life, independence, comfort and mobility of patients.

It allows those affected to get out of bed without help from others – so the independence is maintained and having to leave the familiar surroundings is avoided.

---

## Relief of the nursing staff – more time for the patients

---

**The ROTOFLEX® saves costs and at the same time provides more time to meet the individual needs of patients.**

---

To mobilise those affected and to implement other nursing measures is with the **ROTOFLEX®** not only easier, but also gentler – both for the user and for the nursing staff: transfer and mobilisation of the patient with the stand-up bed is without stress

and easy on the back. Time constraints and physical strain on the nursing staff are significantly reduced. Many measures can be done with less time and personnel costs, thus providing additional time for patient-related activity.

## Getting out of bed – a barrier is overcome

---

**The ROTOFLEX® paves the way back to independence.**

---

Getting out of bed means for people with a movement impairment the overcoming of an invisible barrier. Up till now, they were only able to overcome this with the help of others. The **ROTOFLEX®** overcomes the hurdle, which forces the affected person to lie down in bed and obstructs them from access to activities outside the bed. Four strong, reliable motors lift them up at the touch of a button, program-controlled, softly and quietly for sitting on the edge of the bed and

then back to lying down. All movements of the **ROTOFLEX®** can also be individually and comfortably controlled via the remote control so that the five-part reclining surface can be steplessly brought to any desired position. Nursing staff are relieved considerably by the **ROTOFLEX®**, among other things, with the patient transfer into a wheelchair or a standing aid, body-distant lifting is omitted, down-times are minimised and a longer stay in the occupation becomes possible.



The ROTOFLEX® fulfils all the functions of a fully-fledged nursing bed. It is free-standing, stable and tilt stable. When used in a double bed, the caregiver has access to the patient from all sides, because the lying surface in the turned position can be placed completely horizontally. The stand-up bed is optionally available for left or right rotation.

„Since I had the **ROTOFLEX®** bed, I did not once have to wake my wife in the night. It is a relief for us both. For me, it is a restored sense of life, that I myself can cope alone and do not have to ask for help all the time. My daily nap is also no longer a problem. I go to bed myself and also get out of bed on my own. So, my wife can leave the house at any time during the day.“ *(Bonnich Peterson)*



„The new freedom to organise myself without help, had a positive effect on my battered psyche. In addition, you save money. Before I got the ROTOFLEX®, I needed help from carers at least four times a day. Today the only help I have, is a housekeeper, who does one hour of cleaning every two weeks.“ *(Jan Ronhoff)*

„I can go to bed whenever I want, get up whenever I want, and I do not have nursing staff here around the house all the time. I can do my chores on my own. For me, that means freedom.“ *(Helle Pedersen)*



Patents registered.

GigaFlex®, RotoBett®, RotoCare®, RotoFlex®, PhysioFlex® and PhysioNova® are registered trademarks of PhysioNova GmbH.



# Customisation

As required, the **ROTOFLEX®** is available in different versions up to special construction – the specific requirements for equipment and optics can be realised. It can be

supplemented with additional accessories such as bed lattices, uprights, reversible pressure mattresses and the like and can be expanded as desired and adapted

to the specific requirements of the patient. A detailed overview of the available accessories can be found on our homepage under: [www.physionova.de](http://www.physionova.de).

## Different versions of the ROTOFLEX®



### **ROTOFLEX® P AS A STANDALONE NURSING BED AND ROTOFLEX® E IN THE HOME DOUBLE BED**

The **ROTOFLEX® P** is available as a standalone nursing bed and can also be integrated as an installation bed **ROTOFLEX® E** into the single or double bed at home instead of the existing slatted frame. Visually, the **ROTOFLEX® E** is barely different to a commercial bed.



The **ROTOFLEX® P 250** and the **ROTOFLEX® P 300** are wider than the standard variant. The additionally installed motor gives the reinforced construction more power. The **ROTOFLEX® P 250** and the **ROTOFLEX® P 300** are available between one metre and 1.20 metres in width and create enough space for people with a body weight between 200 and 250 kilogrammes.



The **ROTOFLEX® LOW** is also suitable for small people due to its low-lying surface height. In the lowest position, the lying surface is only 32 cm from the floor – so small people can sit comfortably and safely reach the ground with their feet when standing up. The **ROTOFLEX® LOW** is designed to be particularly easy to dismantle and transport.

The **ROTOFLEX®** has been used in hospitals, retirement homes, rehabilitation centres as well as private homes for more than 20 years. The outstanding functionality of the **ROTOFLEX®** sets new standards in the area of nursing beds, facilitates staff-intensive care procedures and thereby helps to ensure the required quality of care.

### **Do you have questions or would you like to buy a product?**

We can be reached under the phone number **+49 (0) 911 / 977 248 – 0**. We are open from Monday to Friday between 8am and 5pm; our fax number is **+49 (0) 911 / 977 248 – 11**. On our website at [www.physionova.de](http://www.physionova.de) you can also get information about products, accessories and the background of PhysioNova GmbH outside the business hours. Alternatively, you can contact us using the contact form of the homepage or by email to [info@physionova.de](mailto:info@physionova.de).



PhysioNova GmbH  
Vacher Straße 43  
D-90587 Obermichelbach

Tel: +49 (0) 911 / 977 248 – 0  
Fax: +49 (0) 911 / 977 248 – 11  
[info@physionova.de](mailto:info@physionova.de)  
[www.physionova.de](http://www.physionova.de)

