Am Eckfeld II, Gewerbering 8 83543 Rott am Inn

Germany

Tel.: +49 (0)8035-9580448 Fax: +49 (0)8035-9580446

info@hernik.de • www.hernik.de



invoice / delivery address	□ Order
Reference:	Offer

Rehab Car Seat /PA/MAXX

Exportprices Release: 2025

BASIC

Details:

	IPAI-MAX WITHOUT	IPAI-MAX WITH
	extra cushion	extra cushion
SEAT WIDTH (pelvis)	32 cm	26 cm
SEAT WIDTH (knee)	38 cm	38 cm
SEAT DEPTH	33 cm	30 cm
BACK HEIGHT	41 - 56 cm	38 - 53 cm
CROWN HEIGHT	59 - 76 cm	56 - 73 cm
HEIGHT LATERALS	24 - 36 cm	21 - 33 cm
DISTANCE LATERALS	15 - 22 cm	
angle	80° - 120°	
total height	71 - 88 cm	
total width	44 cm	
total depth	44 cm	
weight	from 8,0 kg	

common information (standard):

- crash tested according to ECE-R129/03
- 5-point-harness with upholstery, height adjustable
- belt mechanics installed into the front of the seat
- solid aluminium board, including the mechanism to recline the backrest stepless
- head rest height adjustable, including belt quides
- additional upholstery (20 mm) for the back & headrest
- ISOFIX adapter or clamp bow integrated in base board

IPAI-MAX red-black (450000-E1)

IPAI-MAX grey-black (450000-E2)






accessory that makes the base board NECESSARY


BASE BOARD Necessary for the options below; incl. ISOFIX adapter or clamp bow **ISOFIX** (450000-BEI) (450000-BEK) angular adjustment (5° - 20°) swivel base Conversion set (450100-E) (450200-E) O right from ISOFIX adapter to clamp bow (450020-E) O left from clamp bow to ISOFIX adapter (450030-E) foot rest foot roll (450450-E) swing up* outrigger* (450455-E) (450453-E) (450452-E) *450450-E necessary! *450450-E necessary! further options security pad for the buckle protection mat (to protect the original seats of the car) with snap fastener with magnetic lock (450780-E) (450650-E) (450655-E)




extra cushions



fittings for stabilization





Notes