



1605

## ACCESSORIES FOR SPECIALIST VISION AIDS

# Reading Stand

- The reading stand encourages an upright, orthopaedically correct posture and therefore prevents back pain.
- When the reading stand is used in conjunction with a magnifying vision aid (e.g. system varioPLUS), an optimum angle of view and a uniform distance to the reading material are always achieved, as well as bright illumination of the field of view.
- Guarantees stability and a tremble-free image.
- Supporting surface area: 40 x 35 cm.
- Material: maple wood.
- 5 different angles of adjustment.
- Edge stopper prevents slipping.
- Supplied with clamp attachment for large papers.

### Accessories:

#### 16051 Guide rail

- Elderly people in particular often have problems guiding magnifiers across the reading material. The guide rail facilitates line orientation
- The guide rail can be secured at any height by means of a fixing screw
- Material: birch wood



1605 + 16051