


## Chiroform®™ Sitting Wedge De Luxe

The Chiroform Sitting Wedge De Luxe is made of selected mixture of quality polyurethane chip foam, moulded to 100 KGs pr. Cbm. Combined with heavy duty visco-elastic foam of 65 KGs. Pr. Cbm. The wedge bring you back to the natural position, so that your back suffers no unnecessary strain. In addition, the wedge supports the pelvis and thighs to prevent numbness in the legs. The unique pressure relieving visco-elastic foam has proved its efficiency in the treatment and prevention of decubitus ulcers and back complaints.



### Product numbers

C107/1-B Sitting Wedge De Luxe Black

350x350x65/20 mm

C107/1-G Sitting Wedge De Luxe Grey

350x350x65/20 mm

### Specifications

- Bottom-Foam: 45 to 0 mm bonded chip foam.  
Density (core) ISO 845 / DIN 53420: - approx. 100 Kgs pr. Cbm.
- Contents: Polyetherpolyol, Toluendiisocyanate, Liquid Paraffin, Water, Mine-catalyst, Stannooctoact, Polyethersiloxane.
- Cover: 100% velvet 340 g/m<sup>2</sup>.
- Cleaning: Cover can be hand washed. While washing please turn the cover inside out.

## Practical Advise

Always make sure, that the visco-elastic foam is on the top. Pressing with hand can recognise its slow recovery feeling. The unique visco-elastic foam reacts to body temperature. If room temperature is very low, the visco-elastic foam will become stiff. This is quite normal for the material, which has been constructed to react sensitively on body heat. You will recognise the softening of the material while you "melt" it into your personal shape.

A slightly scent of new material may occur just after removal of transportation packaging, but it should not cause any alarm, as it is harmless and will disappear in the course of a few days.

## Copyrights/Trademarks/Patents

This specification sheet, including, but not limited to, all text, images, figure and the name of Chiroform presented hereon, as well as their arrangement and assembly (collectively, the "Content"), is the property of Chiroform Ltd. and is protected by Danish and international copyright laws. While Chiroform Ltd. makes the Content available to all users, Chiroform Ltd. retains all copyright and related proprietary rights therein. Accordingly, you may not, without the express written permission of Chiroform Ltd.: copy or distribute the Content or any portion thereof; mirror or include the Content or any portion thereof on your own documents; or modify or re-use any portion of the Content. You may make a single copy of the Content for your own personal use. Chiroform Ltd. reserves all other rights in the Content.

Trademark patent no: VA2007 01362 and IR00963673.

## CE Mark

The wedges of Chiroform ®™ conform and complies with the European medical CE mark.

Copyright ®. This Document is the Property of Chiroform, and may not be copied or given to a third part without written permission of Chiroform.



Tingvej 7, DK- 8800 Viborg, Tel. +4586613611  
WWW.CHIROFORM.DK