

# CENTRE FOR TESTING AND CERTIFICATION - MECH-TEST

# **Mechanical Laboratory**

05-077 Warszawa-Wesoła, ul. Klonowa 22 tel.: +48 603 23-26-45, e-mail: cbc.mech.test@gmail.com, www.cbc.org.pl

Date 17.01.2023

# TEST REPORT NO. *CBC* -008/2023

Page 1 of 6

Subject of testing:

Walking aids with built-in handgrips and three or

more legs of which two or more are having wheels.

which provide support whilst walking

Classification according to

Number of specimens: 1

PN-EN ISO 9999:2017-02:12 06 06

Type / Model:

Puma Air rollator

REF: 312501

SN: 0001

Manufacturer:

MOBILEX A/S

Grønlandsvei 5

DK - 8660 Skanderborg

Applicant:

A-Net s.c.

93-469 Łódź.

ul. Łaskowice174

Kind of testing

Mechanical testing for conformity with EN ISO 11199-2 : 2021

Test started: 12.01.2023

Test finished: 17.01.2023

Approved by:

inż. Andrzej Tkaczyk

Special comments / enclosures:

Copyright © 2012 by Centre for Testing and Certification (applicable to report form)

Test results refer only to tested units.

Test results reported here are not applicable to the further modifications of the product affecting its structure, material or technology. This test report shall be neither copied differently as in the whole nor be published without written consent of the Laboratory.



Page: 2 of 6

### CHARACTERISTIC OF PRODUCT

Name: Puma Air rollator Dimension of rollator:

REF: 312501 SN: 0001

Maximum permissible user mass: 9,785 kg 150 kg Mass of rollator: PHOTO OF PRODUCT

	Descripti	on	
E	lements/parameters/materials/d	limensions	Comments
	Distance between	476 mm -	
Dimensions od walking rollator (fig. 2 PN-EN ISO 11199-2)	handgrips (dimension 2)	480 mm	
	Angle between of handgrip	00	
	axis and direction of		
8 =	movement (α)		1
- [토] - 1	Height of rollator	860/818 mm	min.
wa	(dimension 6)	990 mm	max.
요현	Width of rollator	630 mm	
Siz	(dimension 5)		
sio	Turning width	924 mm	
fig	(dimension 1)		
<del> </del>	Length of rollator	764 mm	
-	(dimension 4)		
Dimen	sions of folded rollator (mm)	764 x 850 x 25	56
	Handgrip - diameter	35 mm	Anatomical
ni m			handgrip
Fig.	Handgrip - length	121 mm	
	Front wheels-quantity	2	castor
ato <u>i</u>	Front wheels - diameter	293 mm	wheels
170	Front wheels – width	32/46 mm	
7	Front wheels - brake	none	
S O	Rear wheels - quantity	2	
ee	Rear wheels - diameter	293 mm	
Wheels of rollator	Rear wheels - width	32/46 mm	
	Rear wheels - brake	Included	
Tip	Diameter	none	
	Material	none	
	Colour	none	
	Front legs	Aluminum,	
lla	Bracing member (no. 8)	Hard plastic,	
fr	Rear legs	Bolts, nuts	
10 C	Height adjusting device (no.		
oria fig.	4)		
Material of rollator (fig. 1)	Handgrip (no 5),	Hard plastic	
5 5	Brake elements	_	















≤150 kg 65 cm →I ≨5 kg

CE MD

Page: 3 of 6

Requireme nts according to clause	ne- ding use	Checked	Real	Test result	Comments
Required according to clarity	Test method according to clause	characteristics/assemblies/parameters	Re	Te	Comments
6.1	ISO 14971 ISO 12100	Risk analysis		N/T	
6.2	V/I	Rollators that can be dismantled		N/A	
6.3	V/I	Fasteners		N/A	
6.4	V/I	User mass / load limit	Conf.	Pos.	
6.5	V/I	Equipment of the rollator used inside the premises	Conf.	Pos.	
		• front wheel diameter ≥75mm	Conf.	Pos.	293 mm
		• the presence of two parking brakes	Conf.	Pos.	
6.5	V/I	Equipment of the rollator used outdoors	Conf.	Pos.	
		• front wheel diameter ≥180mm	Conf.	Pos.	293 mm
		• front wheel width ≥22mm	Conf.	Pos.	32 mm
		• the presence of two driving brakes	Conf.	Pos.	
		the presence of two parking brakes	Conf.	Pos.	
6.6		Brakes			
	6.6.1 V/I	Ease of use of the brakes when moving a rollator with more than 2 wheels	Conf.	Pos.	
	6.6.1 V/I	The occurrence of parking brakes in a rollator with more than 2 wheels and a seat for rest, or intended for outdoor use	Conf.	Pos.	
	6.6.1 Measur.	Distance to the brake lever (fig. 3)	Conf.	Pos.	74mm (≤ 75 mm)
	6.6.1 V/I	Ability to compensate for brake wear	Conf.	Pos.	
	6.6.1 V/I	No negative effect on the brake performance when folding, unfolding or adjusting the rollator	Conf.	Pos.	
	6.6.1 V/I	Possibility of adjusting settings (e.g. height) without using tools	Conf.	Pos.	
6.6.2		Brake effectiveness			
	6.6.2.1 Measur.	Parking brake activation/release force	50/30N Conf.	Pos.	≤ 60 N pushing ≤ 40 N pulling
	6.6.2.2	The effectiveness of the travel brake	Conf.	Pos.	Movement of rollator ≤10 mm in 1 mimute
	6.6.2.2	The effectiveness of the parking brake	Conf.	Pos.	Movement of rollator ≤ 10 mm in 1 minute
6.6.3.1	6.6.3.2	Durability of brakes	Conf.	Pos.	100 000 cycles f≤0,5Hz
6.7	Measur.	Handgrip  Recycling	35mm Conf.	Pos.	Width of handgrip ≥20mn and ≤50mm
7.2.1	7.1 ISO 8191-2	Flammability		N/T	
7.2.1	ISO 8191-2			11-	
7.3	ISO 10993-1	Upholstered parts		N/T	
7.4	130 10333-1	Biocompatibility and toxicity		N/T	
7.4	7.4.2	Infection and microbiological contamination		N/T	
7.5	ISO 9227	Cleaning and disinfection  Resistance to corrosion	<u> </u>	N/T	
8	V/I	Ingress of liquids	<del>-</del>	N/T	
9	V/I,			N/A	t <sup>0</sup> ≤ 41 <sup>0</sup> C
,	Weasur.	Temperatures of parts that come in contact with human skin			requirement does not concern heat of direct sold radiation - PN-EN 12182, clause 10a requirement concerns on persons with insensitivenes of skin (who do not feel hed - PN-EN 12182, clause 10
10	10.1, V/I, Measur	Safety of moving parts	Conf.	Pos.	Z. Z
11	11.1, V/I, Measur	Holes and clearances	Conf.	Pos.	
	11.2, V/I, Measur	V-shape openings	Conf.	Pos.	



Mech	anical	Labor	atory	of CBC
IAICOI	antai	1./2117171	ZERUII V	171 1 121

Page: 4 of 6

Requirements according to clause	Test me- thod according to clause	Checked characteristics/assemblies/parameters		Real	Test result	Comments
12	V/I	Folding, adjusting and locking mechanisms		Conf.	Pos.	
	12.3	Locking mechanisms		Conf.	Pos.	
13	V/I, Measur	Carrying handles		·		
		Presence of the handling devices (e.g. hand	les) in		N/A	
		components of mass greater than 10 kg, or			14/24	
		Information indicating the points where cor be lifted and describing how they shall be h disassembly, lifting, carrying and assembly	andled during	-	N/A	
14	V/I	Surfaces, corners and edges		Conf.	Pos.	
15.1	15.2.1	Forward-direction stability		18,8° Conf.	Pos.	Stability required ≥ 15°
15.1	15.2.2	Backward-direction stability		13,9 <sup>0</sup> Conf.	Pos.	Stability required ≥ 7°
15.1	15.2.3	Sideway-direction stability		4,7 ° Conf.	Pos.	Stability required ≥ 3,5
15.1	15.2.4	Stability – with loaded basket,	orwards	17,3° Conf.	Pos.	
		bag, drip, oxygen cylinder	oackwards	17,3 <sup>0</sup> Conf.	Pos.	
			ide	6,0° Conf.	Pos.	
16.1.2	16.1.3	Static strength of resting seat		Conf.	Pos.	1 minute under load 1,2 x user`s weight±2% (180kg)
16.2.2	16.2.3	Static strength of the rollator		Conf.	Pos.	loading 180kg±2%, 1min. NOTE
16.3.2	16.3.3	Strenght of backrest		-	N/A	1,2 x user`s weight±2% x 60sek, x 10 cycles
17.1	17.2	Durability test		Conf.	Pos.	200 000 cycles with load 120kg±2%, f=1Hz
		Information marke		t		
20.2	V/I	a) manufacturer's model identification nam	e and/or numb.	_	N/T	
		b) whether or not the rollator is designed fo outdoor use	r indoor or	-	N/T	
		c) maximum user mass		_	N/T	
		d) name or trade name and address of the m	anufacturer or			
		authorized representative according to local		-	N/T	
		e) year and month of manufacture	1	_	N/T	
		f) maximum safe working load (to be mark accessories	ed on the	_	N/T	
		g) maximum width of the rollator			N/T	
		h) maximum allowed angle between the lor centre line of the handle and the direction of		_	N/T	
		handles are sideways adjustable		_	14/1	
		i) all information shall as far as possible be Pictogram in accordance with ISO 7000 and		_	N/T	
		Instruction	manual			
20.3	V/I	a) information on how to obtain the user into format appropriate for use by people with vor cognitive disabilities		_	N/T	
		b) a description of the intended use including user and the intended environment	ng intended	_	N/T	
		c) maximum user mass		_	N/T	
		d) maximum safe working load for load can		_	N/T	
		accessories such as basket, tray, shopping b				
		e) minimum and maximum height of the ro	llator	_	N/T	
		f) maintenance instructions, if applicable			N/T	
		g) if the rollator is intended to be cleaned, a the method and suitable cleaning materials, precautions needed to avoid corrosion, if ap	including	-	N/T	
		h) if the rollator is intended to be disinfecte of the method and suitable materials, include	d, a description	_	N/T	



Page: 5 of 6

ts according to clause	Test method according to clause	Checked characteristics/assemblies/parameters	Real value	Test result	Comments
		i the overall dimensions (width, length and height) of the rollator, expressed in milimeters, and its mass, expressed in kilograms, when it is ready for use and, if applicable, when it is folded of dismantled	_	N/T	
		j the mass expressed in kilograms, if the rollator can be dismantled or has any removable parts that has a mass that is heavier than 10 kilograms	-	N/T	
		k if the rollator is supposed to be used in combination with other products, the manufacturer shall state which products, and how this can be done in a safe way	_	N/T	
		l a list of accessories, detachable parts and materials that the manufacturer has determined as being intended for use with the rollator	_	N/T	
		m whether and how the rollator can be folded or dismantled to assist in storage or transport	-	N/T	
		n the location and the type of identification number/word on the rollator shall be given for the unique identification number of the assistive product	-	N/T	
		o any adjustment or settings required before the rollator can be used and information on how adjustments or settings affect the rollator	-	N/T	
		p information on adjustment possibilities and the competence required to carry out these adjustments	-	N/T	
		q instructions on operation of all controls	_	N/T	
		r instructions on dismantling and re-assembly of the rollator or any removable parts	_	N/T	
		s the positions of points where the component parts can be gripped for safe moving and handling and/or a method for handling during dismantling, assembly or carrying		N/T	
		t a warning if surface temperatures can increase / decrease when exposed to external sources of heat (e.g. sunlight, outdoor environment	-	N/T	
		u if the intended purpose of the rollator cannot be met without a hazard (e.g. holes, V-shaped opening), a warning and instructions on how to operate the assistive product safety	-	N/T	
		v if the intended purpose of the rollator cannot be met without a hazard due to moving parts such as squeezing, a warning and instructions on how to operate the rollator safety	_	N/T	
		w how to obtain information about the warranty		N/T	
		x warning of the risk of falling from the rollator such as "Incorrect use can lead to hazardous situation - Do not use the products to transport a person".	-	N/T	



Mechanical Laboratory of CBC	Report no.: C	Report no.: CBC-008/2023 Page : 6 of 6			
TEST CONDITIONS					
Ambient temperature	16°C	Required temperature 21°C ±5°C			
Relative humidity of air:	55 %	Not required			
Comments:	*/				
All tests performed with maximum height adjustment of rollator.					
All tests performend in the least stabble position of self-adjusting wheels.					
Tests performed with handles positioned at their maximum (allowed) angl	e to the direction of motion (when adjustmen	t is possible).			
Sequence of tests: stability test, static loading test, fatigue test.					
One rollator was tested.					
During visual inspection before testing any visible defects that cou	ld have influence on test results were not	t stated.			

Pos. – positive; Neg – negative; N/T – not tested; N/A – not applicable; N/R – not required, N/O – not occurred, V/I. – visual inspection, Conf.- conformed.

NOTE 1: Deformation – 24mm, elastic deformation – 22mm, permanent deformation – 2mm (0,21% -  $\leq$ 1%).

### **CONCLUSIONS:**

Testing object conforms with requirements of PN-EN ISO 11199-2: 2021 in scope of mechanical testing ordered by client, excluding biocompatibility tests of material with human body according to PN-EN ISO 10993-1:2010.

---- END ----

