


Following the market requirements and modern design on the health care market we invented the new Leo bed. It meets the highest world's standards and provides full support for disable patients and people who need long term care.

We use a unique pantograph construction to give a flexibility and modern look to the bed. All the electrical items like actuators and control box are hidden under the mattress platform and are not visible for the user. The new bed's construction is independent from the wooden frame, which could be replaced for the new one, without interfering in bed's functions.



Exceptional for Leo bed is enormous height range from 290 mm to 800 mm (Leo 290 version) and 260 mm to 770 mm (Leo 260 version - home care use only) which combines typical, standard bed models with their low variation.

Wide range of adjustment helps medical staff to take a good care of their patients and provides additional safety for the bed users. Solid bed frame gives you a real stability and at the same time a light design.



Leo bed offers many additional equipment for better care giving.

New wooden split side rails are divided 40/60 along the bed. They protect the patient and make lifting and getting up from the bed much easier. They are made of solid wood and silver metal frame which gives feeling of strength and modern design.

Two options of height regulation give opportunity to use the bed at home environment – Leo 290 and also at the health care institutions – Leo 260 for home care use only.

Leo bed offers wide range of additional equipment for example:

- castors with central brake
- detachable foot /head ends
- various shape of the wooden ends
- different kinds of side rails: split side rails, full length side rails, metal side rails
- Battery Back-up System

Modern pantograph construction of bed with 6 points of support gives firm stability. Exceptional for Leo is enormous height range of 50 cm.

Additional abilities are castors with central brake or from two sides individual or in pairs braking option. Castor's central braking option provides wide range of bed's use due to client's needs.

Thanks to this specifications Leo bed combines advantages of beds with lower construction called low with standard version of the bed.



SIDE RAILS - TYPE S

The newest solution – Split Side Rail System. New wooden split side rails are divided 40/60 and protect the patient up to 400 mm of height. According to patients needs there is possibility of lifting part of the side rail from the head side or from the foot side. S Side rails type is a solution which protects the patient during the operation and provides comfort and feeling of patient's safety for medical staff. New Split side rails are manufactured according to current norm EN-60601-2-52.



SIDE RAILS – TYPE L

Full length wooden side rails are the standard equipment of all nursing care beds. They protect the patient up to 385 mm over the mattress platform. For additional protection there is available a set of Side Rails Height Extension, which increases the height of side rails up to 530 mm. This solution is very helpful when the use of pressure mattresses is needed. Full-length wooden side rails are manufactured according to current norm EN-60601-2-52.



SIDE RAILS – TYPE M

The Metal Side Rails system is designed to protect the patient in 3/4th of the bed's length and up to 385 mm of height. Side rails locking system in upper position is operated automatically. Side rails dropping system is not available for the patient and can be only operated by the nurse. Metal Side Rails System is easy and fast operated and doesn't require additional space. Metal side rails are manufactured according to current norm EN-60601-2-52.

BED ENDS

Wide range of bed ends and colours gives opportunity to create your individual style in interior design.

