# **RECUMBENT STEPPER PRO**

# **GB7008**

With half-linked arms and pedals, for focused rehabilitation, user can perform upper and lower body synchronism or independent movements. At the same time, it will strengthen your core muscle groups and enhance your exercise durability.



# **USER ADVANTAGES**























# **CHOOSE YOUR CONSOLE**



Smart Console (10.1" TFT LCD)




Standard Console (9" LCD)













#### A. Patented Independent System

Each independent arm has accommodated three movement pattern including unilateral arms move independently to create single or dual-arm exercises.

# **B.** Direct Wheelchair Access

To meet the specific needs of wheelchair users we have removable seats for wheelchair access.

# **C.** Fully Removable for Wheelchair Access

Simply slide the seat off from the rail for wheelchair access. The removable seat with transportation wheels can make it easily for storage.

#### **D.** Low Step-Through Design

For increasing accessibility.

#### E. Safety Locking for Easy Access

A manual locking mechanism keeps the handle securely in position when you are getting into the seat.

#### F. Multi-Position Ergo Handlebars

Comfortable position for exercising your upper body.