

Pure Air Posture

Clean Air

Completely sealed with clinically proven filter allowing clean air flow, but blocks viruses, allergens, microbes and mold



Soft Touch

Durable, soft touch and high performance textile

Multi Patient Use

Fully waterproof cushions with clean airflow filter. Simply clean and reuse

Self Molding Technology

Microbead technology allows the cushion to adapt to all morphologies forming the patient's body shape

Certified

Certified as Class 1 in EU
Certified as Class 1 within FDA
Approved and registered with the SFDA

Pressure Offloading

Fine microbeads cleverly relieve pressure and promote neutral body alignment


How It Works?

The cushions are composed of high-quality round polystyrene micro particles. Due to the microbead technology the cushion adapts to all morphologies forming the patient's body shape, relieving pressure points. Giving full support and comfort to the patient. Hermetically sealed seams allowing a complete barrier with no stitch holes



PneumaPure Filter Technology Proven Infection Control

Filter allows clean air to flow into the cushions, blocking out pathogens and allergens such as MRSA, C. difficile, E. coli, Pseudomonas aeruginosa, Candida albicans, Aspergillus niger, Norovirus, Influenza Type A, Adenovirus 5 and Coronavirus 229E



Versatile Positioning

Suitable for various applications, including:

- Bed sore prevention and treatment
- Lying positioning during night time
- Preventing asymmetric or poor postures
- Wound care
- Pressure care
- Aligned positioning



The PosiHeel cushion is designed specifically for heel offload. With different internal microbeads best suited for the application.

Designed to be the perfect solution for heel pressure offload.

PosiHeel Cushion

Product No.	Product Name	Qty per Box
PP8054	PosiHeel Offload Positioner 50 x 70cm	6



PosiPerle Cushion Range

Product No.	Product Name	Qty per Box	Image
PP8011	PosiPerle Lateral Support Microbead Cushion 180 x 30cm	3	
PP8012	PosiPerle Half-Moon Microbead Cushion 180dia x 30cm	4	
PP8013	PosiPerle Delta Cushion 60 x 30cm	5	
PP8014	PosiPerle Semi-Fowler Cushion 70 x 56 x 14cm	4	
PP8015	PosiPerle Cylinder Cushion, large 70 x 17cm	12	
PP8016	PosiPerle Universal Cushion, large 60 x 40cm	10	
PP8017	PosiPerle Universal Cushion, small 40 x 30cm	25	
PP8018	PosiPerle Mini Universal Cushion 35 x 22cm	62	
PP8033	PosiPerle Cylinder Cushion, small	25	
PP8034	PosiPerle Abduction Wedge Cushion	25	
PP8035	PosiPerle 1/2 Circle Cushion	16	
PP8045	PosiPerle Semi-Fowler Cushion, small	10	
PP8046	PosiPerle Universal Cushion 45 x 30cm	25	
PP8047	PosiPerle Abduction Wedge Cushion, XS	55	
PP8048	PosiPerle Arc Cushion 90x35cm	12	
PP8049	PosiPerle Arc Cushion 78x30cm	12	
PP8050	PosiPerle Delta Cushion, small 42x25cm	6	
PP8051	PosiPerle Leg Wedge Cushion	4	
PP8055	PosiPerle Small Decubitus 140 x 26cm	3	
PP8056	PosiPerle Small Half-Moon 145 x 24cm	4	
PP8057	PosiPerle Support Wedge 20 x 45 x 15/7.5cm	19	
PP8058	PosiPerle Wedge Cushion 7 x 30 x 20cm	27	
PP8059	PosiPerle XL Delta 40x40x80cm	4	
PP0085	PosiPerle Transport Bag 100 x 30 x 65cm	50	

Covers for PosiPerle Cushions

Cotton and bamboo covers are available for the full PosiPerle range. These are made with an easy access zip opening. Allowing the cushion to be taken out and the cover to be machine washed.

Bamboo Covers

The bamboo covers are made of a more natural wicking, breathable material.



Bamboo Covers Shown

Cotton Covers


The covers are made with an easy access zip opening. Allowing the cushion to be taken out and the cover to be machine washed.



Cotton Cover Shown

How Can I Use PosiPerle?


Head Support



Chest & Side Lying Support



Arm & Hand Support



Leg & Knee Support



Foot & Heel

