

FEMTO S

1.186

USER PROFILE

- ✓ For permanent and long-term fittings
- ✓ Very well suited for a wide range of medical conditions, including the severely afflicted
- ✓ For users with some remaining walking and standing ability (sit to stand transfer)
- ✓ For neurological disorders

THE LIGHT-WEIGHT FOLDING WHEELCHAIR WITH SWING-AWAY AND DETACHABLE LEGRESTS

The FEMTO S is ideal for users looking for a flexible wheelchair suitable for everyday use. The legrests can also be operated without use of the hands, making the FEMTO S particularly suitable for use in institutions or for neurological conditions.


- Active adjustment options for ideal handling
- Ultralight from 9.7 kg
- Frame length depends on the seat depth for optimum ergonomics
- Seat widths from 320 – 520 mm
- Seat depths from 340 – 500 mm
- Optimum adjustability of the castor fork guarantees perfect directional stability
- Easy handling and transfer
- Modular frame concept with hydroformed profile tubes
- Very stable thanks to few screw connections and minimal weight
- Easy adjustment



Reduced frame design



Castor wheel adjustment



Multifunctional legrests




Swing-away and detachable legrests

Folding

Detachable wheels



User weight up to 130 kg



Low empty weight from 9.7 kg up



Seat widths from 320 – 520 mm

